

AIA-360 Assay Specifications

ST FT3 Test Code 076

No.	Item	Data
1	Calib. Req	Show
2	Cal 1	0 pg/mL
3	Cal 2	1.5 pg/mL (example)
4	Cal 3	3.0 pg/mL (example)
5	Cal 4	6.0 pg/mL (example)
6	Cal 5	12 pg/mL (example)
7	Cal 6	25 pg/mL (example)
8	Cal Lot L	
9	Cal Lot R	
10	Unit	pg/mL
11	Smpl. Vol	50
12	Dil. Vol	100
13	Assay L	0.5
14	Assay H	25
15	Ref. L	
16	Ref. H	
17	Decimal	2

Visible only in Test Mode

18	Test Name	#FT3
19	Calib. No	6
20	Calib. Mul	3
21	Calib. Equ	6
22	Calib. CV	90
23	Assay Prtl	1
24	Factor1 A	2.640000
25	Factor1 B	0.000000
26	Factor2 A	1.000000
27	Factor2 B	0.000000
28	V. Conc	0.000000
29	G Origin	0.000000

AIA-900 Assay Specifications

ST FT3 Test Code 076

No.	Item	Data
1	Code	76
2	ACT	1
3	Analyte	#FT3
4	Lot	***Enter current cal lot no.
5	CAL 1	0 pg/mL
6	CAL 2	1.5 pg/mL (example)
7	CAL 3	3.0 pg/mL (example)
8	CAL 4	6.0 pg/mL (example)
9	CAL 5	12 pg/mL (example)
10	CAL 6	25 pg/mL (example)
11	Cal lot L	
12	Cal lot R	
13	Unit	pg/mL
14	Decimal	2
15	Assay Low	0.5
16	Assay High	25
17	Reference Low	
18	Reference High	
19	Reschedule Low	0.5
20	Reschedule High	25
21	Factor A	1
22	Factor B	0
23	Sample Volume	50
24	Diluent Volume	100
25	2 Step reagent dispensing volume	0
26	Calibration code	45
27	CAL. No.	6
28	CAL. MUL.	3
29	CAL. EQU.	6
30	CAL. CV	90
31	DIL. SP1	1
32	DIL. SP2	1
33	DIL. CAL. (Calculation of dil ratio of conc.)	1
34	DIL. CNTL.	1
35	DIL. DO	1
36	DIL. AH.	1
37	DIL. CALC.	1
38	DIL. CODE	0
39	DIL. NAME	
40	DIL. PRTY	3
41	PRE. SPVOL (Vol. of pretreated sample)	0
42	PRE. 1VOL (Vol. of pretreatment sol-1)	0
43	PRE. 2 VOL (Vol. of pretreatment sol-2)	0
44	PRE. CODE (pretreated sol. code)	0
45	PRE. NAME (pretreated sol. name)	0
46	Protocol	1
47	SYS. F_A	2.64
48	SYS. F_B	0
49	V. CONC.	0
50	G. ORIGIN	0

AIA-2000 Assay Specifications

ST FT3 Test Code 076

No.	Item	Data
1	Unit	pg/mL
2	Decimal places	2
3	Reference low	
4	Reference high	
5	Reschedule low	0.5
6	Reschedule high	25
7	Assay range low	0.5
8	Assay range high	25
9	Specimen diluent code	
10	Specimen diluent name	
11	Dilution factor for Sp. 1	1
12	Dilution factor for Sp. 2	1
13	Dilution factor for Control	1
14	Default multiplier for DO	1
15	Default multiplier for >H	1
16	Factor A	1.00000e+000
17	Factor B	0.00000e+000
18	Calibration code	6
19	Calibrator replicates	3
20	Calibrator lot	***Enter current cal. lot no
21	Calibrator Concentration	
22	Cal - 01	0 pg/mL
23	Cal - 02	1.5 pg/mL (example)
24	Cal - 03	3.0 pg/mL (example)
25	Cal - 04	6.0 pg/mL (example)
26	Cal - 05	12 pg/mL (example)
27	Cal - 06	25 pg/mL (example)
28		
29		
30		
31		
32		
33		
34	Dilution factor for calibrator	1
35	Calibration period	90
36	Incubation time (10/40)	10
37	Assay protocol	1
38	Specimen volume	50
39	Diluent volume	100
40	Conjugate volume	0
41	Diluted conjugate volume	0
42	Pretreatment	0
43	Pretreatment specimen volume	0
44	Pretreatment Reagent code	0
45	Pretreatment Reagent name	0
46	Pretreatment Reagent 1 volume	0
47	Pretreatment Reagent 2 volume	0
48	System Factor	2.64
49	Calibration Code Check	045
50	Virtual Concentration	0.0
51	Graph Origin	0.0
52	Calculation with Dilution Factor	Yes

Free Triiodothyronine ST AIA-PACK FT3

Name and Intended Use

ST AIA-PACK FT3 is designed for IN VITRO DIAGNOSTIC USE ONLY for the quantitative measurement of Free Triiodothyronine (FT3) in human serum or heparinized plasma on specific Tosoh AIA System analyzers.

Summary and Explanation of Test

Triiodothyronine (T3) is present in human serum in an equilibrium mixture of bound and free forms, with approximately 0.4% of the total T3 circulating as Free Triiodothyronine (FT3). Any change in the serum concentration of binding proteins will cause a parallel rise in the concentration of total T3 with the FT3 remaining relatively unchanged.¹ Its physiological action is apparent only when free fractions of T3 and T4 are taken in by target cell membrane receptors.²⁻⁴ Since protein bound T3 is not available for cellular uptake and does not appear to exert any metabolic effect, the measurement of FT3 may therefore represent a more valid clinical assessment of patient thyroid status.⁵⁻¹⁰

Direct measurement of free T3 enables thyroid function examination even in the presence of abnormal liver function, hormone fluctuation during pregnancy,¹¹ and variations in levels of serum binding proteins.

Principle of the Assay

The ST AIA-PACK FT3 is a competitive enzyme immunoassay which is performed entirely within the AIA-PACK. Free Triiodothyronine (FT3) present in the test sample competes with enzyme-labeled triiodothyronine (T3) for a limited number of binding sites on a T3-specific antibody immobilized on magnetic beads. The beads are washed to remove the unbound enzyme labeled free triiodothyronine and are then incubated with a fluorogenic substrate, 4-methylumbelliferyl phosphate (4MUP). The amount of enzyme labeled free triiodothyronine that binds to the beads is inversely proportional to the free triiodothyronine concentration in the test sample. A standard curve using a range of known standard concentrations is constructed and unknown free triiodothyronine concentrations are calculated using this curve.

Material Provided (ST AIA-PACK FT3, Cat. No. 0025286)

5 trays x 20 test cups (ST AIA-PACK FT3 Test Cup)

Plastic test cups containing lyophilized magnetic beads with anti-T3 sheep monoclonal antibody and T3 conjugated to bovine alkaline phosphatase with 0.1% sodium azide as a preservative.

Materials Required But Not Provided

The following materials are not provided but are required to perform Free Triiodothyronine analysis using the ST AIA-PACK FT3 (Cat. No. 0025286) on specific Tosoh AIA Systems. They are available separately from Tosoh.

Materials	Cat. No.
AIA-SYSTEMS:	
AIA-360	0019945
AIA-900	0022930
AIA-900 9tray Sorter	0022931
AIA-900 19tray Sorter	0022932
AIA-2000 ST	0022100
AIA-2000 LA	0022101
AIA-PACK:	
AIA-PACK Substrate Set II	0020968
AIA-PACK Substrate/Reconstituent	
AIA-PACK FT3 Calibrator Set	0020386
Calibrator #1 0 pg/dL	
Calibrator #2 1.5 pg/dL (approx.)	
Calibrator #3 3.0 pg/dL (approx.)	
Calibrator #4 6.0 pg/dL (approx.)	
Calibrator #5 12 pg/dL (approx.)	
Calibrator #6 25 pg/dL (approx.)	
AIA-PACK Wash Concentrate Set	0020955
AIA-PACK Diluent Concentrate Set	0020956
AIA-PACK Detector Standardization Test Cups	0020970
Sample Cups	0018581
ADDITIONAL REQUIREMENTS: (Except AIA-360)	
Pipette Tips (1000/Pkg)	019215
Tip Rack (Empty)	019216
Preloaded Pipette Tips (96 Tips X 50 Racks)	996010
Preloaded Pipette Tips (96 Tips X 5 Racks)	996005

Only materials obtained from Tosoh should be used. Materials obtained elsewhere should not be substituted since assay performance is based strictly on Tosoh materials.

Warnings and Precautions

- The ST AIA-PACK FT3 is intended for in vitro diagnostic use only.
- **Rx** only. US Federal law restricts this device to sale by or on the order of a licensed healthcare practitioner.
- Inspect the packaging and the exterior of the aluminum pouch or the vials for any sign of damage before use. If any damages are visible, contact your local Tosoh sales representative.
- Test cups from different lots should not be mixed within a tray.
- The ST AIA-PACK FT3 contains sodium azide, which may react with lead or copper plumbing to form potentially explosive metal azides. When disposing of such reagents, always flush with large volumes of water to prevent azide build-up.
- Human sera is not used in the preparation of this product, however, since human specimens will be used for samples and other quality control products in the lab may be derived from human serum, please use standard laboratory safety procedures in handling all specimens and controls.
- For safe waste disposal, it is recommended that each laboratory complies with established laboratory procedures and local, state, and federal regulations.
- Do not use beyond the expiration date.
- Serum, dust, metal, or microorganism contamination may cause degradation of reconstituted substrate solution. Store in a clean environment, away from direct sunlight and ultraviolet light.
- Tosoh Automated Immunoassays utilizing alkaline phosphatase-based technologies should not be used with samples from patients under Asfotase Alfa treatment.

Storage and Stability

All unopened materials are stable until the expiration date on the label when stored at the specified temperature.

Materials	Cat. No.
Refrigerator Temperature (2 - 8 °C):	
ST AIA-PACK FT3	0025286
AIA-PACK FT3 Calibrator Set	0020386
AIA-PACK Substrate Set II	0020968
AIA-PACK Wash Concentrate	0020955
AIA-PACK Diluent Concentrate	0020956
Room Temperature (18 - 25 °C):	
AIA-PACK Detector Standardization Test Cups	0020970

ST AIA-PACK FT3 test cups may be stored for up to 24 hours at room temperature of 18 - 25 °C. Calibrators must be kept tightly sealed and refrigerated at 2 - 8 °C. After opening, calibrators should be used within 24 hours. Reconstituted substrate solution is stable for 3 days at 18 - 25 °C

or 30 days at 2 - 8 °C. Working diluent and wash solutions are stable for 30 days at room temperature (18 - 25 °C). Reagents should not be used if they appear cloudy or discolored.

Specimen Collection and Handling

Serum or heparinized plasma is required for the assay. EDTA and citrated plasma **SHOULD NOT BE USED**.

No special patient preparation is necessary. When using serum, a venous blood sample is collected aseptically without additives. Store at 18 - 25 °C until a clot has formed (usually 15 - 45 minutes), then centrifuge to obtain the serum specimen for assay.

To use heparinized plasma, a venous blood sample is collected aseptically with the designated additive. Centrifuge and separate plasma from the packed cells as soon as possible.

Samples may be stored at 2 - 8 °C for up to 24 hours prior to analysis. If the analysis cannot be done within 24 hours, the sample should be stored frozen at -20 °C or below for up to 60 days.

Repeated freeze-thaw cycles should be avoided. Turbid serum samples or samples containing particulate matter should be centrifuged prior to testing. Prior to assay, bring frozen samples to room temperature (18 - 25 °C) slowly and mix gently.

The sample required for analysis is 50 µL.

Procedure

1) Reagent Preparation

1a) Substrate Solution

Bring all reagents to room temperature (18 - 25 °C) before preparing the working reagent. Add the entire contents of the Substrate Reconstituent (100 mL) to the lyophilized Substrate and mix thoroughly to dissolve the solid material.

1b) Wash Solution

Add the entire contents of the Wash Concentrate (100 mL) to approximately 2.0 L of CAP Class I water or the clinical laboratory reagent water, mix well, and adjust the final volume to 2.5 L.

1c) Diluent

Add the entire contents of the Diluent Concentrate (100 mL) to approximately 4.0 L of CAP Class I water or the clinical laboratory reagent water, mix well, and adjust the final volume to 5.0 L.

2) Calibration


2a) Calibration Curve

The calibrators for use with the ST AIA-PACK FT3 are prepared gravimetrically and are compared to internal reference standards.

The calibration curve for the ST AIA-PACK FT3 is stable for up to 90 days. Calibration stability is monitored by quality control performance and is dependent on proper reagent handling and AIA System maintenance according to the manufacturer's instructions.

Recalibration may be necessary more frequently if controls are out of the established range for this assay or if certain service procedures are performed (e.g. temperature adjustment, sampling mechanism changes, or detector lamp adjustment or change). For further information regarding instrument operation, consult the AIA Systems Operator's Manual.

The sample calibration curve below shows the algorithm used to calculate the results. This is an example only. Actual results will vary depending on the Instrument type and lot number used.



2b) Calibration Procedure

Refer to the appropriate AIA System Operator's Manual for procedural instructions.

- Verify that both the calibrator lot and concentration numbers have been correctly entered into the software.
- The AIA-PACK FT3 CALIBRATOR SET is lyophilized. All levels should be reconstituted with 1.0 mL of CAP Class I water or the clinical laboratory reagent water.
- Tosoh recommends that all calibrators be run in triplicate.

2c) Calibration Acceptability criteria

- i) Since there is an inverse relationship between concentration and rate, the rates should decrease as the concentration increases.
- ii) The replicate values should be within a 10% range.

2d) Calibration Review and Acceptance

- i) Review the calibration curve carefully, using the criteria listed above.
- ii) Edit the calibration if necessary, then accept the calibration.
- iii) For further information regarding calibration, consult the specific AIA System Operator's Manual.

3) Quality Control

3a) Commercially Available Controls

Commercially available controls should be run at least once per day. It is recommended that at least two (2) levels of controls, normal and abnormal, be used. Laboratory policy for this particular assay designates the following:

Control Material: _____
Frequency: _____

Lot number of control material, acceptable limits, and corrective action to be taken if controls do not meet laboratory criteria will be found in a separate quality control document maintained by the laboratory.

3b) Quality Control Procedure

- i) Assay quality control specimens as instructed in the specific Operator's Manual for your analyzer. In addition, refer to the AIA System Operator's Manual for detailed instructions on defining and editing the files.
- ii) Quality control material to be run with this assay is defined by individual laboratory policy.

4) Specimen Processing

4a) Preparation

Following specific instructions in the Operator's Manual for the analyzer, place samples on the instrument appropriately.

4b) Assay Procedure

- i) Ensure a sufficient quantity of ST AIA-PACK FT3 test cups for the number of samples to be run.
- ii) Load patient samples as instructed in the Operator's Manual and proceed with analysis.

Procedural Notes

1. Lyophilized Substrate must be completely dissolved.
2. Ligand assays performed by the Tosoh AIA System Analyzer require that the laboratory use water designated by the CAP Class I or by the clinical laboratory reagent water. Water should be tested at least once per month and should be free of particulate matter including bacteria. The pH of the water should also be routinely tested. For further information, consult the CLSI document GP40-A4-AMD, Preparation and Testing of Reagent Water in the Clinical Laboratory; Approved Guideline - Fourth Edition.
3. If a specimen free triiodothyronine concentration is found to be greater than 25 pg/mL, results should be reported as greater than 25 pg/mL. Dilution or recovery studies cannot be performed on the Free Triiodothyronine assay. The equilibrium between bound and free hormone in the Free Triiodothyronine assay is altered both by dilutions and the addition of measured amounts of T3.
4. The AIA systems can store two different calibration curves for each analyte at one time. Therefore, up to two different lots of ST AIA-PACK FT3 Test cups can be used during the same run.
5. If the assay specifications for this test are not already in the system software, the specifications must be entered under test code 076.

Calculation of Results

The AIA Systems perform all sample and reagent handling operations automatically. The AIA Systems read the rate of fluorescence produced by the reaction and automatically convert the rate to Free Triiodothyronine concentration in pg/mL.

Evaluation of Results

Quality Control

In order to monitor and evaluate the precision of the analytical performance, it is recommended that commercially available control samples be assayed daily.

The minimum recommendations for the frequency of running internal control material are:

- After calibration, two levels of controls are run in order to accept the calibration curve.
- The two levels of controls are also repeated after calibration when certain service procedures are performed (e.g. temperature adjustment, sampling mechanism changes, maintenance of the wash probe or detector lamp adjustment or change).
- After daily maintenance, at least two levels of the control should be run in order to verify the overall performance of the Tosoh AIA System Analyzers.

If one or more control sample value(s) is out of the acceptable range, it will be necessary to investigate the validity of the calibration curve before reporting patient results.

Standard laboratory procedures should be followed in accordance with the regulatory agency under which the laboratory operates.

Limitations of the Procedure

For diagnostic purposes, the results obtained from this assay should be used in conjunction with other data (e.g., symptoms, results of other tests, clinical impressions, therapy, etc.).

Using ST AIA-PACK FT3, the highest concentration of Free Triiodothyronine measurable is approximately 25.0 pg/mL, and the lowest measurable concentration is 0.5 pg/mL (assay sensitivity).

Although the approximate value of the highest calibrator is 25.0 pg/mL, the exact concentration may be slightly different. The assay specification, Assay Range High, should be defined as the upper limit of the assay range, 25 pg/mL.

Although hemolysis has an insignificant effect on the assay, hemolyzed samples may indicate mistreatment of a specimen prior to assay and results should be interpreted with caution.

Lipemia has an insignificant effect on the assay except in the case of gross lipemia where spatial interference may occur.

Heparin may cause in vivo effects in free T3 assays. Blood collection for free T3 assays should not be performed concurrent with administration of heparin therapy.

Certain medications may interfere with assay performance. Specimens from patients taking medicines and/or medical treatment may show erroneous results. All results should be interpreted with respect to the clinical picture of the patient.

Specimens from patients under treatment with Diclofenac sodium and Mefenamic acid may demonstrate falsely elevated results. Patients with anti-T3 antibody may also show false results.

Samples from patients who had an injection of fluorescein, which is used in fluorescein fundus angiography, may cause falsely elevated results.

For a more complete understanding of the limitations of this procedure, please refer to the Specimen Collection and Handling, Warnings and Precautions, Storage and Stability, and Procedural Notes sections in this insert sheet.

Expected Values

Each laboratory should determine a reference interval which corresponds to the characteristics of the population being tested. As with all diagnostic procedures, clinical results must be interpreted with regard to concomitant medications administered to the patient.¹²

Reference Ranges


The following data was collected using the AIA-PACK FT3. The ST AIA-PACK FT3 demonstrates equivalent performance. (See Comparative Analysis.)

Two studies were performed using apparently healthy individuals.

In the first, a reference interval representing the central 95% of values tested with the AIA-PACK FT3 was determined in serum samples from 820 apparently healthy individuals. Mean value was 3.1 pg/mL.

Reference interval: 2.0 - 4.9 pg/mL (3.1 - 7.5 pmol/L)

In the second study, 49 apparently healthy individuals were tested for FT3 using the AIA-PACK FT3 assay. Mean value was 2.83 pg/mL with a standard deviation of 0.35 pg/mL.



Conversion Factors

FT3 concentrations in this application are in units of pg/mL. Conversion to SI units of pmol/L may be made using the following equation:


$$\text{pmol FT3/L} = \text{pg FT3/mL} \times 1.536$$

Performance Characteristics

The following performance characteristics were determined using the Tosoh AIA NexIA Automated Immunoassay Analyzer. The AIA-600 II demonstrates equivalent performance.

1) Comparative Analysis

Samples were collected at random. Results from the ST AIA-PACK FT3 assay (y) were compared to those obtained using the AIA-PACK FT3 assay (x). A total of 44 serum samples were tested in the study.



Comparative Analysis
ST AIA-PACK FT3 vs AIA-PACK FT3

Slope:	1.0378
y-intercept:	-0.2156
Correlation Coefficient:	0.9964
Range of Samples:	0.83 – 28.89 pg/mL
n:	44

2) Precision

2a) Intra-assay precision

The intra-assay (within run) precision coefficient of variation was determined using three controls in a total of 20 runs. Within each run, one set of duplicates per control was assayed. The mean of each duplicate was used to obtain the pooled standard deviation (SD), which was then used to calculate the coefficient of variation (CV).

Sample	Mean (pg/mL)	Pooled SD (pg/mL)	Coefficient of Variation (%)
Serum A	2.648	0.124	4.7
Serum B	5.439	0.212	3.9
Serum C	13.572	0.321	2.4

2b) Inter-assay precision

Inter-assay precision was evaluated at three different concentrations by analyzing samples in 20 different runs.

Sample	Number of Replicates	Mean (pg/mL)	Standard Deviation (pg/mL)	Coefficient of Variation (%)
Serum A	20	2.716	0.149	5.5
Serum B	20	5.471	0.240	4.4
Serum C	20	13.708	0.544	4.0

2c) Total precision

Total precision was determined by the duplicate assay of three controls in 20 separate runs. The means of each run were used to calculate the standard deviation (SD) and coefficient of variation (CV).

Sample	Mean (pg/mL)	Pooled SD (pg/mL)	Coefficient of Variation (%)
Serum A	2.648	0.183	6.9
Serum B	5.439	0.262	4.8
Serum C	13.572	0.492	3.6

Specificity

The following substances were tested for cross-reactivity. The cross-reactivity (%) is the percent of the compound which will be identified as FT3. If these compounds are present in the specimen at the same concentration as FT3, the final result will be increased by these percentages.

Compound	Amount Added (pg/mL)	Amount Measured (pg/mL)	Cross-reactivity (%)
L-Triiodothyronine (L-T3)	15.00	15.00	100
D-Triiodothyronine (D-T3)	15.00	18.75	125
L-Thyroxine (L-T4)	15.00	0.0225	0.15
D-Thyroxine (D-T4)	15.00	0.0105	0.07
3, 3', 5'-Triiodothyronine	15.00	0.0450	0.30
3, 5-Diiodo-L-tyrosine	15.00	0.000	<0.001

Sensitivity

The minimal detectable concentration (MDC) of Free Triiodothyronine is estimated to be 0.5 pg/mL. The MDC is defined as that concentration of FT3 which corresponds to the rate of fluorescence that is two standard deviations from the mean rate of fluorescence of 20 replicate determinations of a zero calibrator.

Interference

Interference is defined, for purposes of this study, to be recovery outside of 10% of the known specimen mean concentration.

- Added hemoglobin (up to 390 mg/dL), free bilirubin (up to 17 mg/dL) and conjugated bilirubin (up to 19 mg/dL) do not interfere with the assay.
- Lipemia, as indicated by added triglyceride (up to 830 mg/dL), does not interfere with the assay.
- Ascorbic acid (up to 20 mg/dL) does not interfere with the assay.
- Tosoh Automated Immunoassays utilizing alkaline phosphatase-based technologies should not be used with samples from patients under Asfotase Alfa treatment.

References

1. Burtis, C.A. and Ashwood, E.R. [edited by], 1994, Tietz Textbook of Clinical Chemistry, 2nd Edition, W.B Saunders, 1718.
2. Oppenheimer, J.H., 1986, Role of Plasma Proteins in the Binding, Distribution and Metabolism of the Thyroid Hormones. New Eng. J. Med., 278:1153.
3. Prince, M.P. and Ramsden, D.B., 1977, A New Theoretical Description of the Binding of Thyroid Hormones by Serum Proteins. Clin. Endocrinol., 7:307.
4. Chopra, I.J., et al, 1978, Pathway of Metabolism of Thyroid Hormones. Rec. Prog. Horm. Res., 34: 521.
5. Ekins, R.P., 1984, Free Hormones in Blood: The Concept and the Measurement. J. Clin. Immunoassay, 7:163.
6. Sterling, K. and Lazarus, J.H., 1977, The Thyroid and its Control. Ann. Rev. Physiol., 39: 349.
7. Borst, G.C., et al, 1983, Euthyroid Hyperthyroxinemia. Ann. Int. Med., 98:366.
8. Utiger, R.D., 1974, Serum Triiodothyronine in Man. Ann. Rev. Med., 25: 289.
9. Lerman, J., 1953, Physiologic Activity of L-Triiodothyronine. J. Clin. Endocrinol., 13:1341.
10. Nauman, J.A., Nauman, A., Werner, S.C., 1967, Total and Free Triiodothyronine in Human Serum. J. Clin. Invest., 46:1346.
11. Sparre, L.S., et al, 1987, Oestrogen and Thyroxine Binding Globulin Levels in Early Normal Pregnancy. Acta Endocrinol., 114: 298.
12. Young, D., 1990, Effects of Drugs on Clinical Laboratory Tests. 3rd Edition, Washington, DC, American Association for Clinical Chemistry Press.



TOSOH CORPORATION
Bioscience Division
 2-2-1 Yaesu, Chuo-ku, Tokyo 104-0028 Japan
 Phone: +81 3 6636 3734
 Fax: +81 3 6636 3627
 Website: www.tosoh.com

Supplied by

TOSOH BIOSCIENCE, INC.
 3600 Gantz Road,
 Grove City, OH 43123 USA
 Phone: +1 650 615 4970
 Fax: +1 650 615 0415
 Phone: (800) 248 6764
 Fax: (800) 685 7595
 Website: www.tosohbioscience.us


European Conformity


In vitro diagnostic medical device


Consult instructions for use


Temperature limitation


Batch code / Lot number


Manufacturer


Authorized representative
in the European Community


Use by date


Catalogue number
/ Part number


Supplied by


Sufficient for

US: For prescription use only

AIA-PACK

FT3 Calibrator Set

Intended Use

The AIA-PACK FT3 Calibrator Set is intended for IN VITRO DIAGNOSTIC USE ONLY for the calibration the ST AIA-PACK FT3 Assay.

Summary and Explanation

The AIA-PACK FT3 Calibrator Set contains human serum with assigned levels of free triiodothyronine (FT3). Calibration should be performed according to the schedule indicated in the AIA System Operator's Manual. The calibrators in this set are prepared gravimetrically and are compared to internal reference standards.

Materials Provided (Cat. No. 0020386)

2 x 1 mL	Calibrator	No. 1	0 pg/mL	
			Human serum containing no detectable concentration of FT3 (0 pg FT3/mL), and 0.1% sodium azide as a preservative. (Lyophilized)	
2 x 1 mL	Calibrator	No. 2	1.5 pg/mL	(approx.)
		No. 3	3.0 pg/mL	(approx.)
		No. 4	6.0 pg/mL	(approx.)
		No. 5	12.0 pg/mL	(approx.)
		No. 6	25.0 pg/mL	(approx.)
			Human serum containing the assigned concentration of FT3 (described on each vial), and 0.1% sodium azide as a preservative. (Lyophilized)	

Warnings and Precautions

- The AIA-PACK FT3 Calibrator Set is for in vitro diagnostic use.
- These materials contain sodium azide, which may react with lead or copper plumbing to form potentially explosive metal azides. When disposing of such reagents, always flush with large volumes of water to prevent azide build-up.
- Human sera used in the preparation of this product has been tested by FDA cleared methods and found negative for the presence of HBsAg and antibody to HIV-1 and HCV. Because no test method can offer complete assurance that products derived from human blood will not transmit infectious agents, it is recommended that this product be handled with the same precautions as used for patient samples.
- Do not use beyond the expiration date.

Preparation and Storage

- Using a volumetric pipette or a calibrated adjustable pipette, reconstitute the lyophilized calibrators accurately with 1.0 mL of CAP Class I water or the clinical laboratory reagent water. Allow the lyophilized material to fully dissolve, then mix the calibrators gently but thoroughly prior to performing the calibration.
- Bring calibrator to room temperature (18 - 25 °C) for use.
- Always store the Calibrator Set in an upright position at 2 - 8 °C when not in use.

Stability

When stored unopened and refrigerated at 2 - 8 °C, the AIA-PACK FT3 Calibrator Set is stable until the expiration date on the label. After opening, the calibrators should be used within 24 hours.

Procedure

Refer to the CALIBRATION PROCEDURE in the AIA-PACK section of this analyte application. For additional procedural instructions regarding calibration, refer to the AIA System Operator's Manual.

1. Verify that both the calibrator lot and concentrations have been correctly entered into the analyzer software.
2. Add the appropriate amounts of each calibrator to sample cups (refer to the instrument worklist for the amount required in each sample cup).
3. Place the sample cups and the test cups on the analyzer as indicated on the worklist.
4. Tosoh recommends that all calibrators be run in triplicate.

Assignment of Values

The AIA-PACK FT3 Calibrator Set contains assigned concentrations of Free Triiodothyronine. The assigned value is determined on a lot-to-lot basis and is designed to provide an assay calibration range of approximately 0 to 25 pg/mL of Free Triiodothyronine. The Calibrator Set has been prepared gravimetrically and is compared to internal reference standards.

Results

- Since there is an inverse relationship between concentration and rate, the rates should decrease as the concentration increases.
- The replicate values should be within a 10% range.

Limitations

The ST AIA-PACK FT3 Calibrator Set is designed solely for use with AIA-PACK assay procedures.

References

1. AIA Analyte Application Manual (AAM). Tosoh Bioscience, Inc., Grove City, OH.
2. AIA-System Operator's Manual. Tosoh Bioscience, Inc., Grove City, OH.



TOSOH CORPORATION
Bioscience Division
2-2-1 Yaesu, Chuo-ku, Tokyo 104-0028 Japan
Phone: +81 3 6636 3734
Fax: +81 3 6636 3627
Website: www.tosoh.com

Supplied by

TOSOH BIOSCIENCE, INC.
3600 Gantz Road,
Grove City, OH 43123 USA
Phone: +1 650 615 4970
Fax: +1 650 615 0415
Phone: (800) 248 6764
Fax: (800) 685 7595
Website: www.tosohbioscience.us


European Conformity


In vitro diagnostic medical device


Consult instructions for use


Temperature limitation


Batch code / Lot number


Manufacturer


Authorized representative
in the European Community


Use by date


Catalogue number
/ Part number

Supplied by




Sufficient for

