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SE-HPLC; size-exclusion HPLC, Non-red CE; non-reduced capillary electrophoresis

Ref. Y. Wen et al., Asian J. Complementary and Alternative Medicine, 2025, 4, <http://dx.doi.org/10.53043/2347-3894.acam13002> を一部改変

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- 東ソー、テクニカルノート (トヨパール) No. 21、「二重特異性抗体の精製は、これで決まり！」
- Y. Li et al., A roadmap for IgG-like bispecific antibody purification, Approaches to the Purification, Analysis and Characterization of Antibody-Based Therapeutics, 2020, p.167-179, <https://doi.org/10.1016/B978-0-08-103019-6.00008-4>
- S. W. Chen et al., Current trends and challenges in the downstream purification of bispecific antibodies, Antibody Therapeutics, 2021, 4, 73-88, <https://doi.org/10.1093/abt/tbab007>
- Y. Li, General strategies for IgG-like bispecific antibody purification, Biotechnol. Prog. 41 (2025) 41:e3515, <https://doi.org/10.1002/btpr.3515>
- M. S. Karbyshev et al., Trends and challenges in bispecific antibody production, J. Chromatogr., A, 1744 (2025) 465722, <https://doi.org/10.1016/j.chroma.2025.465722>
- Q. Li et al., The downstream purification of bispecific antibodies, Anal. Biochem., 696 (2025) 115692, <https://doi.org/10.1016/j.ab.2024.115692>
- K. Mix et al., Rapid affinity-based purification of multi-specific antibodies using Kappa Select and Protein L, MABS 17 (2025) 2483272, <https://doi.org/10.1080/19420862.2025.2483272>
- N. Ingavat et al., Affinity resin selection for efficient capture of bispecific antibodies as guided by domain composition, Process Biochemistry 154 (2025) 1-11, <https://doi.org/10.1016/j.procbio.2025.04.007>
- Y. Wen et al., Removal of Antibody-Related Fragments During Asymmetric Bispecific Antibody Purification, Asian Journal of Complementary and Alternative Medicine, 2025, 4, <http://dx.doi.org/10.53043/2347-3894.acam13002>

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