

AIA-360 Assay Specifications

ST AFP Test Code 014

No.	Item	Data
1	Calib. Req	Show
2	Cal 1	0 ng/mL
3	Cal 2	200 ng/mL (example)
4		
5		
6		
7		
8	Cal Lot L	
9	Cal Lot R	
10	Unit	ng/mL
11	Smpl. Vol	25
12	Dil. Vol	100
13	Assay L	1.0
14	Assay H	400.0
15	Ref. L	
16	Ref. H	
17	Decimal	1

Visible only in Test Mode

18	Test Name	#AFP
19	Calib. No	2
20	Calib. Mul	3
21	Calib. Equ	1
22	Calib. CV	90
23	Assay Prtl	1
24	Factor1 A	1.000000
25	Factor1 B	0.000000
26	Factor2 A	1.000000
27	Factor2 B	0.000000
28	V. Conc	0.000000
29	G Origin	0.000000

AIA-900 Assay Specifications

ST AFP Test Code 014

No.	Item	Data
1	Code	14
2	ACT	1
3	Analyte	#AFP
4	Lot	***Enter current cal. lot no.
5	CAL 1	0 ng/mL
6	CAL 2	200 ng/mL (example)
7		
8		
9		
10		
11	Cal lot L	
12	Cal lot R	
13	Unit	ng/mL
14	Decimal	1
15	Assay Low	1.0
16	Assay High	400.0
17	Reference Low	
18	Reference High	
19	Reschedule Low	1.0
20	Reschedule High	400.0
21	Factor A	1
22	Factor B	0
23	Sample Volume	25
24	Diluent Volume	100
25	2 Step reagent dispensing volume	0
26	Calibration code	2
27	CAL. No.	2
28	CAL. MUL.	3
29	CAL. EQU.	1
30	CAL. CV	90
31	DIL. SP1	1
32	DIL. SP2	1
33	DIL. CAL. (Calculation of dil ratio of conc.)	1
34	DIL. CNTL.	1
35	DIL. DO	10
36	DIL. AH.	5
37	DIL. CALC.	1
38	DIL. CODE	2
39	DIL. NAME	AFP
40	DIL. PRTY	3
41	PRE. SPVOL (Vol. of pretreated sample)	0
42	PRE. 1VOL (Vol. of pretreatment sol-1)	0
43	PRE. 2 VOL (Vol. of pretreatment sol-2)	0
44	PRE. CODE (pretreated sol. code)	0
45	PRE. NAME (pretreated sol. name)	0
46	Protocol	1
47	SYS. F_A	1
48	SYS. F_B	0
49	V. CONC.	0
50	G. ORIGIN	0

AIA-2000 Assay Specifications

ST AFP Test Code 014

No.	Item	Data
1	Unit	ng/mL
2	Decimal places	1
3	Reference low	
4	Reference high	
5	Reschedule low	1.0
6	Reschedule high	400.0
7	Assay range low	1.0
8	Assay range high	400.0
9	Specimen diluent code	2
10	Specimen diluent name	AFP
11	Dilution factor for Sp. 1	1
12	Dilution factor for Sp. 2	1
13	Dilution factor for Control	1
14	Default multiplier for DO	100
15	Default multiplier for >H	10
16	Factor A	1.00000e+000
17	Factor B	0.00000e+000
18	Calibration code	1
19	Calibrator replicates	3
20	Calibrator lot	***Enter current cal. lot no.
21	Calibrator Concentration	
22	Cal - 1	0 ng/mL
23	Cal - 2	200 ng/mL (example)
24		
25		
26		
27		
28		
29		
30		
31		
32		
33		
34	Dilution factor for calibrator	1
35	Calibration period	90
36	Incubation time (10/40)	10
37	Assay protocol	1
38	Specimen volume	25
39	Diluent volume	100
40	Conjugate volume	0
41	Diluted conjugate volume	0
42	Pretreatment	0
43	Pretreatment specimen volume	0
44	Pretreatment Reagent code	0
45	Pretreatment Reagent name	0
46	Pretreatment Reagent 1 volume	0
47	Pretreatment Reagent 2 volume	0
48	System Factor	1
49	Calibration Code Check	002
50	Virtual Concentration	0.0
51	Graph Origin	0.0
52	Calculation with Dilution Factor	Yes

Alpha-Fetoprotein ST AIA-PACK AFP

Caution: The sale and distribution of this device is restricted by United States federal law to, by, or on the order of a physician. In addition, the use of this device is restricted to, by, or on the order of a physician. Because of differences in reagent specificity and assay methods, the concentration of AFP in a given specimen may vary with devices from different manufacturers. Values obtained with different assay methods cannot be used interchangeably. It is mandatory that results reported by the laboratory to the physician include the identity of the assay used. If the assay method for AFP is changed during the course of monitoring patients with serial AFP levels, baseline values for the patients being serially monitored must be confirmed by additional sequential testing.

Name and Intended Use

ST AIA-PACK AFP is designed for IN VITRO DIAGNOSTIC USE ONLY for the quantitative measurement of Alpha-Fetoprotein (AFP) in serum and plasma to aid in the management of patients with nonseminomatous testicular carcinoma.

Summary and Explanation of Test

Alpha-fetoprotein (AFP) is a fetal globulin synthesized by the fetal liver, gastrointestinal tract and yolk sac.¹ AFP is a glycoprotein with a molecular weight of 70,000 daltons containing 4.3% carbohydrate, of which approximately 20% is sialic acid.² Within a year after birth, AFP concentration in a normal serum decreases to a barely detectable level.³ Initially the test was used as an aid in the diagnosis of hepatomas,⁴ since increased serum AFP concentrations were first observed in some patients with primary hepatocellular carcinoma. Elevated serum AFP levels may also be observed in patients with nonseminomatous testicular carcinoma.⁵⁻⁷ Increased serum AFP levels have not been observed in patients with testicular carcinoma which is of purely seminoma origin. In both nonseminomatous testicular carcinoma and primary hepatocellular carcinoma in which AFP elevations were observed, the levels frequently decline during successful therapy.⁸⁻¹⁰ Because of this, serial monitoring of serum AFP levels in such patients has been shown to facilitate the clinical management of the disease.

A few other common cancers, usually with hepatic metastases, increase AFP levels; these include carcinoma of stomach, gallbladder, prostate, and bronchi. Increased AFP levels have also been reported in non-malignant diseases including ataxia telangiectasia, hereditary tyrosinemia, neonatal hepatitis, atresia of the biliary tract,¹¹⁻¹⁴ tumor of other tissues¹⁵⁻¹⁷ as well as in pregnant women.^{18,19}

AFP measurements are not recommended for use as a general screening procedure to detect the presence of cancer in the general population. Likewise, since not all tumors are AFP producing, a serum AFP measurement within the reference intervals found in apparently healthy subjects does not rule out the possibility of the presence of malignant disease.

Principle of the Assay

The ST AIA-PACK AFP is a two-site immunoenzymometric assay which is performed entirely in the AIA-PACK. AFP present in the test sample is bound with monoclonal antibody immobilized on a magnetic solid phase and enzyme-labeled monoclonal antibody in the AIA-PACK. The magnetic beads are washed to remove unbound enzyme-labeled monoclonal antibody and are then incubated with a fluorogenic substrate, 4-methylumbelliferyl phosphate (4MUP). The amount of enzyme-labeled monoclonal antibody that binds to the beads is directly proportional to the AFP concentration in the test sample. A standard curve is constructed, and unknown sample concentrations are calculated using this curve.

Material Provided (ST AIA-PACK AFP, Cat. No. 0025252)

5 trays x 20 test cups (ST AIA-PACK AFP Test Cup)

Plastic test cups containing lyophilized magnetic beads coated with anti-AFP mouse monoclonal antibody and mouse monoclonal antibody (to human AFP) conjugated to bovine alkaline phosphatase with 0.1% sodium azide as a preservative.

Materials Required But Not Provided

The following materials are not provided but are required to perform Alpha-Fetoprotein analysis using the ST AIA-PACK AFP (Cat. No. 0025252) on the Tosoh AIA Systems. They are available separately from Tosoh.

	Materials	Cat. No.
AIA-SYSTEMS:		
	AIA-360	0019945
	AIA-900	0022930
	AIA-900 9tray Sorter	0022931
	AIA-900 19tray Sorter	0022932
	AIA-2000 ST	0022100
	AIA-2000 LA	0022101
AIA-PACK:		
	AIA-PACK Substrate Set II	0020968
	AIA-PACK Substrate/Reconstituent	
	AIA-PACK AFP Calibrator Set	0020352
	Calibrator Zero 0 ng/mL	
	Positive 200 ng/mL (approx.)	
	AIA-PACK AFP Sample Diluting Solution	0020552
	AIA-PACK Wash Concentrate Set	0020955
	AIA-PACK Diluent Concentrate Set	0020956
	AIA-PACK Detector Standardization Test Cups	0020970
	AIA-PACK Sample Treatment Cups	0020971
	Sample Cups	0018581
ADDITIONAL REQUIREMENTS: (Except AIA-360)		
	Pipette Tips (1000/Pkg)	019215
	Tip Rack (Empty)	019216

Preloaded Pipette Tips (96 Tips X 50 Racks)
Preloaded Pipette Tips (96 Tips X 5 Racks)

996010
996005

Only materials obtained from Tosoh should be used. Materials obtained elsewhere should not be substituted since assay performance is based strictly on Tosoh materials.

Warnings and Precautions

- The ST AIA-PACK AFP is intended for in vitro diagnostic use only.
- **Rx** only. US Federal law restricts this device to sale by or on the order of a licensed healthcare practitioner.
- Inspect the packaging and the exterior of the aluminum pouch or the vials for any sign of damage before use. If any damages are visible, contact your local Tosoh sales representative. Test cups from different lots should not be mixed within a tray.
- The ST AIA-PACK AFP contains sodium azide, which may react with lead or copper plumbing to form potentially explosive metal azides. When disposing of such reagents, always flush with large volumes of water to prevent azide build-up.
- Human sera is not used in the preparation of this product, however, since human specimens will be used for samples and other quality control products in the lab may be derived from human serum, please use standard laboratory safety procedures in handling all specimens and controls.
- For safe waste disposal, it is recommended that each laboratory complies with established laboratory procedures and local, state, and federal regulations.
- Do not use beyond the expiration date.
- The ST AIA-PACK AFP has been designed so that the high dose “hook effect” is not a problem for the vast majority of samples. Samples with AFP concentrations between 400 and 200,000 ng/mL will read > 400 ng/mL. The “hook effect” phenomenon may occur at AFP concentrations > 200,000 ng/mL.
- Serum, dust, metal, or microorganism contamination may cause degradation of reconstituted substrate solution. Store in a clean environment, away from direct sunlight and ultraviolet light.
- Tosoh Automated Immunoassays utilizing alkaline phosphatase-based technologies should not be used with samples from patients under Asfotase Alfa treatment.

Storage and Stability

All unopened materials are stable until the expiration date on the label when stored at the specified temperature.

Materials	Cat. No.
Refrigerator Temperature (2 - 8 °C):	
ST AIA-PACK AFP	0025252
AIA-PACK AFP Calibrator Set	0020352
AIA-PACK Sample Diluting Solution	0020552

AIA-PACK Calibration Verification/Linearity Test Set	0020652
AIA-PACK Substrate Set II	0020968
AIA-PACK Wash Concentrate	0020955
AIA-PACK Diluent Concentrate	0020956

Room Temperature (18 - 25 °C):

AIA-PACK Detector Standardization Test Cups	0020970
AIA-PACK Sample Treatment Cups	0020971

ST AIA-PACK AFP test cups may be stored for up to 24 hours at a room temperature of 18 - 25 °C. Calibrators must be kept tightly sealed and refrigerated at 2 - 8 °C. After opening, calibrators should be used within 24 hours. After opening, Sample Diluting Solution is stable for up to 90 days refrigerated at 2 - 8 °C. Reconstituted substrate solution is stable for 3 days at 18 - 25 °C or 30 days at 2 - 8 °C. Working diluent and wash solutions are stable for 30 days at room temperature (18 - 25 °C). Reagents should not be used if they appear cloudy or discolored.

Specimen Collection and Handling

Serum or heparinized plasma is required for the assay. EDTA and citrated plasma **SHOULD NOT BE USED**.

No special patient preparation is necessary. When using serum, a venous blood sample is collected aseptically without additives. Store at 18 - 25 °C until a clot has formed (usually 15 - 45 minutes), then centrifuge to obtain the serum specimen for assay.

To use heparinized plasma, a venous blood sample is collected aseptically with the designated additive. Centrifuge and separate plasma from the packed cells as soon as possible.

Samples may be stored at 2 - 8 °C for up to 24 hours prior to analysis. If the analysis cannot be done within 24 hours, the sample should be stored frozen at -20 °C or below for up to 60 days.

Repeated freeze-thaw cycles should be avoided. Turbid serum samples or samples containing particulate matter should be centrifuged prior to testing. Prior to assay, slowly bring frozen samples to room temperature (18 - 25 °C) and mix gently.

The sample required for analysis is 25 µL.

Procedure

1) Reagent Preparation

1a) Substrate Solution

Bring all reagents to room temperature (18 - 25 °C) before preparing the working reagent. Add the entire contents of the Substrate Reconstituent (100 mL) to the lyophilized Substrate and mix thoroughly to dissolve the solid material.

1b) Wash Solution

Add the entire contents of the Wash Concentrate (100 mL) to approximately 2.0 L of CAP Class I water or the clinical laboratory reagent water, mix well, and adjust the final volume to 2.5 L.

1c) Diluent

Add the entire contents of the Diluent Concentrate (100 mL) to approximately 4.0 L of CAP Class I water or the clinical laboratory reagent water, mix well, and adjust the final volume to 5.0 L.

2) Calibration


2a) Calibration Curve

The calibrators for use with the ST AIA-PACK AFP have been standardized against WHO 1st IS 72/225 (1974).

The calibration curve for the ST AIA-PACK AFP is stable for up to 90 days. Calibration stability is monitored by quality control performance and is dependent on proper reagent handling and AIA System maintenance according to the manufacturer's instructions.

Recalibration may be necessary more frequently if controls are out of the established range for this assay or if certain service procedures are performed (e.g. temperature adjustment, sampling mechanism changes, or detector lamp adjustment or change). For further information regarding instrument operation, consult the AIA System Operator's Manual.

The sample calibration curve below shows the algorithm used to calculate the results. This is an example only. Actual results will vary depending on the Instrument type and lot number used.



2b) Calibration Procedure

Refer to the appropriate AIA System Operator's Manual for procedural instructions.

- i) Verify that the calibrator concentrations listed on the vial labels as well as the calibrator lot numbers have been correctly entered into the software.
- ii) Calibrators for ST AIA-PACK AFP are provided ready to use.
- iii) Tosoh recommends that all calibrators be run in triplicate.

2c) Calibration Acceptability criteria

- i) The mean rate for the zero calibrator should be < 3.0 nM/sec.
- ii) Since there is a direct relationship between concentration and rate, the rates should increase as the concentration increases.
- iii) The replicate values should be within a 10% range.

2d) Calibration Review and Acceptance

- i) Using the criteria above, review the calibration curve carefully.
- ii) Edit the calibration if necessary, then accept the calibration.
- iii) For further information regarding calibration, consult the specific AIA System Operator's Manual.

3) Quality Control

3a) Commercially Available Controls

Commercially available controls should be run at least once per day. It is recommended that at least two (2) levels of controls, normal and abnormal, be used. Laboratory policy for this particular assay designates the following:

Control Material: _____
Frequency: _____

Lot number of control material, acceptable limits, and corrective action to be taken if controls do not meet laboratory criteria will be found in a separate quality control document maintained by the laboratory.

3b) Quality Control Procedure

- i) Assay quality control specimens as instructed in the specific Operator's Manual for your analyzer. In addition, refer to the AIA System Operator's Manual for detailed instructions on defining and editing the files.
- ii) Quality control material to be run with this assay is defined by individual laboratory policy.

4) Specimen Processing

4a) Preparation

Following specific instructions in the Operator's Manual for the analyzer, place samples on the instrument appropriately.

4b) Assay Procedure

- i) Ensure a sufficient quantity of ST AIA-PACK AFP test cups for the number of samples to be run.
- ii) Load patient samples as instructed in the Operator's Manual and proceed with analysis. Note: The AIA-900 and AIA-2000 will require AIA-PACK Sample Treatment Cups if onboard dilutions are utilized.

Procedural Notes

1. Lyophilized Substrate must be completely dissolved.
2. Ligand assays performed by the Tosoh AIA System Analyzer require that the laboratory use water designated by the CAP Class I or by the clinical laboratory reagent water. Water should be tested at least once per month and should be free of particulate matter including bacteria. The pH of the water should also be routinely tested. For further information, consult the CLSI document GP40-A4-AMD, Preparation and Testing of Reagent Water in the Clinical Laboratory; Approved Guideline - Fourth Edition.
3. If a serum specimen Alpha-Fetoprotein concentration is found to be greater than the linearity limit of the assay, 400 ng/mL, the specimen should be diluted with the AIA-PACK AFP Sample Diluting Solution and reassayed according to the Assay Procedure. The recommended dilution for samples containing greater than 400 ng/mL is 1:10 or 1:100. It is desirable to dilute the serum sample so that the diluted sample reads between 10 and 400 ng/mL. For further information on the dilution of specimens, refer to the AIA System Operator's Manual.
4. The AIA systems can store two different calibration curves for each analyte at one time. Therefore, up to two different lots of ST AIA-PACK AFP Test cups can be used during the same run.
5. If the assay specifications for this test are not already in the system software, the specifications must be entered under test code 014.

Calculation of Results

The AIA Systems perform all sample and reagent handling operations automatically. The AIA Systems read the rate of fluorescence produced by the reaction and automatically convert the rate to Alpha-Fetoprotein concentration in ng/mL.

For detailed information regarding programming dilutions, consult the appropriate Operator's Manual.

Evaluation of Results

Quality Control

In order to monitor and evaluate the precision of the analytical performance, it is recommended that commercially available control samples be assayed daily.

The minimum recommendations for the frequency of running internal control material are:

- After calibration, at least two levels of controls should be run in order to accept the calibration curve.
- The two levels of controls should also be repeated after calibration when certain service procedures are performed (e.g. temperature adjustment, sampling mechanism changes, maintenance of the wash probe or detector lamp adjustment or change).
- After daily maintenance, at least two levels of the control should be run in order to verify the overall performance of the Tosoh AIA System Analyzers.

If one or more control sample value(s) is out of the acceptable range, it is necessary to investigate the validity of the calibration curve before reporting patient results.

Standard laboratory procedures should be followed in accordance with the regulatory agency under which the laboratory operates.

Limitations of the Procedure

For diagnostic purposes, the results obtained from this assay should be used in conjunction with other data (e.g., symptoms, results of other tests, clinical impressions, therapy, etc.).

Using ST AIA-PACK AFP, the highest concentration of Alpha-Fetoprotein measurable without dilution is approximately 400 ng/mL, and the lowest measurable concentration is 1.0 ng/mL (assay sensitivity).

Although the approximate value of the highest calibrator is 200 ng/mL, the exact concentration may be slightly different. The assay specification, Assay Range High, should be defined as the upper limit of the assay range, 400 ng/mL.

Although hemolysis has an insignificant effect on the assay, hemolyzed samples may indicate mistreatment of a specimen prior to assay and results should be interpreted with caution.

Lipemia has an insignificant effect on the assay except in the case of gross lipemia where spatial interference may occur.

Specimens from patients who have received preparations of mouse monoclonal antibodies for diagnosis or therapy may contain human anti-mouse antibodies (HAMA). Such specimens may show falsely elevated values when tested for Alpha-Fetoprotein.

Certain medications may interfere with assay performance. Specimens from patients taking medicines and/or medical treatment may show erroneous results. All results should be interpreted with respect to the clinical picture of the patient.

For a more complete understanding of the limitations of this procedure, please refer to the Specimen Collection and Handling, Warnings and Precautions, Storage and Stability, and Procedural Notes sections in this insert sheet.

Expected Values

Each laboratory should determine a reference interval which corresponds to the characteristics of the population being tested. As with all diagnostic procedures, clinical results must be interpreted with regard to concomitant medications administered to the patient.²⁰

Reference Ranges

The interval given here was determined in serum samples from 204 domestic adult, non-pregnant, and apparently healthy individuals.

Reference Interval = ≤ 5.63 ng/mL (≤ 4.65 IU/mL)

Category	Samples	ng/mL	IU/mL
Male	104	≤ 6.13 ng/mL	(≤ 5.07 IU/mL)
Female	100	≤ 4.98 ng/mL	(≤ 4.12 IU/mL)

Conversion Factors

AFP concentrations in this application are in units of ng/mL. Conversion to SI units of IU/mL may be made using the following equation:

$$\text{IU AFP/mL} = \text{ng AFP/mL} \div 1.21$$

The following data was collected using the AIA-PACK AFP. The ST AIA-PACK AFP demonstrates equivalent performance. (See the Comparative Analysis.)

The serum concentration of AFP was measured with AIA-PACK AFP in a series of monitoring studies drawn from patients diagnosed with non-seminomatous germ cell carcinoma. The patients were separated into six categories. A brief description of the categories with examples of patient series from each category is presented below.

- Category 1: Initial high AFP. Orchiectomy followed by decreasing AFP. NED with AFP in normal range.
- Category 2: Initial high AFP. AFP decreasing post surgery/chemotherapy to mirror clinical condition.
- Category 3: Initial high AFP history. AFP in normal range post surgery/chemotherapy. Remaining low with NED.
- Category 4: Recurrence of high AFP with metastases. AFP decreasing with chemotherapy and remission.
- Category 5: No history of AFP testing at initial stage of disease. AFP normal with NED.
- Category 6: No elevation of AFP pre-orchiectomy.

Summary of Serial Monitoring with AIA-PACK AFP in Patients with Nonseminomatous Testicular Carcinoma

	Initial Value ng/mL	0-3 months ng/mL	3-6 months ng/mL	6-12 months ng/mL	12-18 months ng/mL	18-24 months ng/mL	24-36 months ng/mL	36-48 months ng/mL	+48 months ng/mL
Category 1									
Patient 1	121	119 138 400 525	57.0 26.0 6.0	3.0 4.0	n/a	6.0	4.0	4.0	4.0
Patient 2	113,000	61,222 3,327 713 26.0	16.0 15.0	8.0 6.0	4.0 5.0	6.0	5.0 5.0	n/a	n/a
Category 2									
Patient 3	15.0	9.0 11.0 92.0	37.0	n/a	n/a	n/a	n/a	n/a	n/a
Patient 4	258	278 378 190 29.0 18.0 11.0	8.0 7.0	n/a	3.0	n/a	n/a	n/a	n/a
Category 3									
Patient 5	n/a	1.0	n/a	1.0	1.0	n/a	1.0	1.0	1.0
Patient 6	n/a	1.0	n/a	2.0	1.0	1.0	1.0	2.0	2.0
Category 4									
Patient 7	97.0	271 783	57.0 9.0 21.0 24.0	16.0 16.0 12.0 11.0	n/a	5.0	4.0	3.0	2.0
Patient 8	2.0	3.0	17.0	4.0 3.0	3.0 2.0	2.0	2.0 2.0	2.0 2.0	2.0
Category 5									
Patient 9	n/a	2.0	n/a	n/a	n/a	1.0	2.0	n/a	2.0
Patient 10	n/a	4.0	3.0	3.0	4.0	n/a	n/a	n/a	n/a
Category 6									
Patient 11	n/a	2.0	2.0 2.0	2.0 2.0	1.0 1.0	n/a	1.0 1.0	1.0	n/a
Patient 12	n/a	2.0	n/a	2.0	2.0 2.0	2.0 1.0	2.0	n/a	n/a

Distribution of Serum AIA-PACK AFP Concentrations in Healthy Subjects, Malignant and Benign Disease States

	Number of Samples	0 - 6.0 ng/mL	6.1- 10.0 ng/mL	10.1- 20.0 ng/mL	20.1- 100.0 ng/mL	100.1- 300.0 ng/mL	> 300 ng/mL
Healthy Subjects							
Men	104	96%	4%	0%	0%	0%	0%
Women	100	98%	2%	0%	0%	0%	0%
Total	204	97%	3%	0%	0%	0%	0%
Carcinoma							
Testicular							
Non-seminoma	50	52%	4%	6%	18%	8%	12%
Seminoma	11	91%	9%	0%	0%	0%	0%
Primary							
Hepatocellular	51	20%	2%	2%	14%	4%	58%
Gastrointestinal	20	90%	10%	0%	0%	0%	0%
Genitourinary	31	81%	16%	0%	3%	0%	0%
Ovarian	20	85%	15%	0%	0%	0%	0%
Pancreatic	5	100%	0%	0%	0%	0%	0%
Other	20	95%	5%	0%	0%	0%	0%
Nonmalignant Diseases							
Cirrhosis	19	89.5%	0%	10.5%	0%	0%	0%
Hepatitis	40	75%	10%	7.5%	2.5%	5%	0%
Mammary	10	90%	0%	10%	0%	0%	0%
Pancreatitis	10	100%	0%	0%	0%	0%	0%
Pulmonary	10	80%	10%	10%	0%	0%	0%
Other	10	90%	10%	0%	0%	0%	0%

In this study, 97% of the healthy subjects had serum AFP concentrations less than or equal to 6.0 ng/mL. Each laboratory should establish its own normal range. Serum AFP concentrations as determined using the ST AIA-PACK AFP assay should be used in conjunction with other clinical and laboratory data.

Performance Characteristics

The following performance characteristics were established using the Tosoh AIA NexIA Automated Immunoassay Analyzer.

1) Accuracy

1a) Recovery:

Three serum pools were spiked with three different levels of AFP and assayed before and after spiking.

Sample	Initial Value (ng/mL)	AFP Added (ng/mL)	Expected Value (ng/mL)	Measured Value (ng/mL)	Percent Recovery (%)
Serum A1	2.63	95.11	97.73	96.40	98.6
	2.63	190.21	192.84	186.47	96.7
	2.63	380.42	383.05	380.20	99.3
Serum B1	4.85	95.11	99.96	92.45	92.5
	4.85	190.21	195.05	193.04	99.0
	4.85	380.42	385.27	375.72	97.5
Serum C1	8.79	95.11	103.89	101.21	97.4
	8.79	190.21	199.00	200.97	101.0
	8.79	380.42	389.21	390.86	100.4

1b) Dilution:

Three serum samples containing high concentrations of AFP were serially diluted with AIA-PACK AFP Sample Diluting Solution and assayed.


Sample	Dilution Factor	Expected Value (ng/mL)	Measured Value (ng/mL)	Percent Recovery (%)
Serum A2	none		353.83	
	7.5/10	265.37	259.45	97.8
	5.0/10	176.92	171.45	96.9
	2.5/10	88.46	82.49	93.3
	1.0/10	35.38	34.64	97.9
Serum B2	none		383.15	
	7.5/10	287.36	285.96	99.5
	5.0/10	191.57	191.22	99.8
	2.5/10	95.79	94.17	98.3
	1.0/10	38.32	38.83	101.3
Serum C2	none		346.35	
	7.5/10	259.76	256.39	98.7
	5.0/10	173.17	172.11	99.4
	2.5/10	86.59	86.07	99.4
	1.0/10	34.63	35.59	102.8

1c) Method Comparison:

Samples were collected at random. Results from the ST AIA-PACK AFP (y) were compared to those obtained using the AIA-PACK AFP (x). A total of 50 serum samples were tested in the study.

Slope 0.972 (95% CI: 0.963 to 0.980)
y-intercept 0.011 (95% CI: -1.2490 to 1.2709)
Correlation Coefficient 0.9995
Range of Samples (x) 1.313 – 414.112
Number of Samples (n) 50

Figure 1.
Comparative Analysis
ST AIA-PACK AFP vs AIA-PACK AFP (40 min)



2) Precision

2a) Intra-assay precision

The intra-assay (within run) precision was determined using three controls in a total of 20 runs. Within each run, one set of duplicates per control was assayed. The mean of each duplicate was used to obtain the pooled standard deviation (SD), which was then used to calculate the coefficient of variation.

Sample	Mean (ng/mL)	Pooled SD (ng/mL)	Coefficient of Variation (%)
Serum A3	10.63	0.24	2.3
Serum B3	102.32	2.17	2.1
Serum C3	286.82	3.85	1.3

2b) Inter-assay precision

The inter-assay (between run) precision coefficient of variation was evaluated at three different concentrations by analyzing samples in 20 separate runs.

Sample	Number of Replicates	Mean (ng/mL)	Standard Deviation (ng/mL)	Coefficient of Variation (%)
Serum A3	20	10.27	0.38	3.7
Serum B3	20	99.67	2.91	2.9
Serum C3	20	279.46	7.71	2.8

2c) Total precision

Total precision was determined by the duplicate assay of three controls in 20 separate runs. The means of each run were used to calculate the standard deviation (SD) and coefficient of variation.

Sample	Mean (ng/mL)	Standard Deviation (ng/mL)	Coefficient of Variation (%)
Serum A3	10.63	0.58	5.5
Serum B3	102.32	4.47	4.4
Serum C3	286.82	12.54	4.4

Sensitivity

The minimal detectable concentration (MDC) of Alpha-Fetoprotein is estimated to be 1.0 ng/mL. The MDC is defined as that concentration of AFP which corresponds to the rate of fluorescence that is two standard deviations from the mean rate of fluorescence of 20 replicate determinations of a zero calibrator.

Interference

Interference is defined, for purposes of this study, to be recovery outside of 10% of the known specimen mean concentration.

- Added hemoglobin (up to 390 mg/dL), free bilirubin (up to 17 mg/dL) and conjugated bilirubin (up to 19 mg/dL) do not interfere with the assay.
- Lipemia, as indicated by added triglyceride (up to 1,660 mg/dL), does not interfere with the assay.
- Protein, as indicated by added albumin (up to 5 g/dL), does not interfere with the assay.
- Ascorbic acid (up to 20 mg/dL) does not interfere with the assay.
- Chemotherapeutic Agents: Therapeutic concentrations of the following chemotherapeutic agents added to serum samples which were then assayed, do not interfere with ST AIA-PACK AFP Assay: Cisplatin, Velban, Bleomycin, Adriamycin, Vincristine, Methotrexate, 5-Fluorouracil, and Mitomycin C.
- Tosoh Automated Immunoassays utilizing alkaline phosphatase-based technologies should not be used with samples from patients under Asfotase Alfa treatment.

References

1. Gitlin, D., et al., 1972, Synthesis of Alpha-Fetoprotein by Liver, Yolk Sac, and Gastrointestinal Tract of the Human Conceptus. *Cancer Res*, 32:979.
2. Ruoslahti, E., Seppala, M., 1971, Studies of Carcino-Fetal Proteins: Physical and Chemical Properties of Human Alpha-Fetoprotein. *Int J Cancer*, 7:218.
3. Masseyeff, R., et al., 1974, Evolution en Fonction de l'Age du Taux Serique Physiologique de l'Alpha-Fetoprotein chez l'Homme et le Rat. (French) in: *Alpha-Feto-Protein* (Masseyeff, R., ed.) p.313. INSERM, Paris.
4. Tatarinov, Y.S., 1964, Finding of an Embryospecific Alpha- Globulin in Blood Stream in Patient with Primary Hepatic Cancer (Russian). *Vopr Med Khim*;10:90.
5. Javadpour, N., et al., 1978, Human Chorionic Gonadotropin (HCG) and Alpha-Fetoprotein (AFP) in Sera and Tumor Cells of Patients with Testicular Seminoma. A Prospective Study. *Cancer*, 42:2768.
6. Lange P.H., et al., 1980, Serum Alpha-Fetoprotein and Human Chorionic Gonadotropin in Patients with Seminoma. *J Urol*,124:472.
7. Jacobsen, G.K., 1983, Alpha-Fetoprotein (AFP) and Human Chorionic Gonadotropin (HCG) in Testicular Germ Cell Tumours. *Acta Path. Microbiol. Immunol Scand Sect A*, 91:183.
8. Anderson, T., et al., 1979, Testicular Germ-Cell Neoplasms: Recent Advances in Diagnosis and Therapy. *Ann Int Med*, 90:373.
9. Bosl, G.J., et al., 1981, Human Chorionic Gonadotropin and Alphafetoprotein in the Staging of Nonseminomatous Testicular Cancer. *Cancer*, 47:328.
10. Vogelzang, N.J., et al., 1982, Acute Changes of Alpha-Fetoprotein and Human Chorionic Gonadotropin during Induction Chemotherapy of Germ Cell Tumors. *Cancer Res*, 42:4855.


11. Waldmann, T.A., McIntire, K.R., 1972, Serum-Alpha-Fetoprotein Levels in Patients with Ataxia-Telangiectasia. *Lancet*, 2:1112.
12. Belanger, L., 1973, Tyrosinémie Héritaire et Alpha-1-Fetoprotéine. II. Recherche Tissulaire Comparée de l'Alpha-Fetoprotéine dans Deux Cas de Tyrosinémie Héritaire. Considérations sur l'Ontogenèse de la Foetoprotéine Humaine (French). *Pathol Biol*, 21:457.
13. Kew, M.C., et al., 1973, Serum Alpha-Fetoprotein Levels in Acute Viral Hepatitis. *Gut*, 14:939.
14. Endo, Y., et al., 1975, Clinical Significance of Alpha-Fetoprotein in Hepatitis and Liver Cirrhosis. *Ann N Y Acad Sci*, 259:234.
15. van Nagell, J.R., et al., 1981, Biochemical Markers in the Plasma and Tumors of Patients with Gynecologic Malignancies. *Cancer*, 48:495.
16. Maeyama, M., et al., 1984, Serial Serum Determination on Alpha-Fetoprotein as a Marker of the Effect of Postoperative Chemotherapy in Ovarian Endodermal Sinus Tumor. *Gynecol Oncol*, 17:104.
17. Yasunami, R., et al., 1981, Primary Lung Cancer Producing Alpha- Fetoprotein: A Case Report. *Cancer*, 47:926.
18. Purves, L.R., Purves, M., 1972, Serum Alpha-Fetoprotein. VI. The Radio-Immunoassay Evidence for the Presence of AFP in the Serum of normal People and During Pregnancy. *S Afr Med J*, 46:1290.
19. Seppala, M., Ruoslahti, E., 1973, Alpha Fetoprotein: Physiology and Pathology During Pregnancy and Application to Antenatal Diagnosis. *J Perinat Med*, 1:104.
20. Young, D., 1990, Effects of Drugs on Clinical Laboratory Tests. 3rd Edition, Washington, DC, American Association for Clinical Chemistry Press.



TOSOH CORPORATION
Bioscience Division
 2-2-1 Yaesu, Chuo-ku, Tokyo 104-0028 Japan
 Phone: +81 3 6636 3734
 Fax: +81 3 6636 3627
 Website: www.tosoh.com

Supplied by

TOSOH BIOSCIENCE, INC.
 3600 Gantz Road,
 Grove City, OH 43123 USA
 Phone: +1 650 615 4970
 Fax: +1 650 615 0415
 Phone: (800) 248 6764
 Fax: (800) 685 7595
 Website: www.tosohbioscience.us

 European Conformity	 In vitro diagnostic medical device	 Consult instructions for use	 Temperature limitation
 Batch code / Lot number	 Manufacturer	 Authorized representative in the European Community	 Use by date
 Catalogue number / Part number	 Supplied by	 Sufficient for	 US: For prescription use only

AIA-PACK

AFP Calibrator Set

Intended Use

The AIA-PACK AFP Calibrator Set is intended for IN VITRO diagnostic use only for the calibration of the ST AIA-PACK AFP Assays.

Summary and Explanation

The AIA-PACK AFP Calibrator Set contains human serum with assigned levels of Alpha-Fetoprotein. Calibration should be performed according to the schedule indicated in the AIA System Operator's Manual. The calibrators in this set have been standardized against WHO 1st IS 72/225 (1974).

Material Provided (Cat. No. 0020352)

2 x 1 mL Zero Calibrator

Human serum containing no detectable concentration of AFP (0 ng AFP/mL), and 0.1% sodium azide as preservative (ready to use).

2 x 1 mL Positive Calibrator

Human serum containing the assigned concentration of AFP (approximately 200 ng/mL, described on each vial) and 0.1% sodium azide as preservative (ready to use).

Warnings and Precautions

- The AIA-PACK AFP Calibrator Set is for in vitro diagnostic use.
- These materials contain sodium azide, which may react with lead or copper plumbing to form potentially explosive metal azides. When disposing of such reagents, always flush with large volumes of water to prevent azide build-up.
- Human sera used in the preparation of this product has been tested by FDA cleared methods and found negative for the presence of HBsAg and antibody to HIV-1 and HCV. Because no test method can offer complete assurance that products derived from human blood will not transmit infectious agents, it is recommended that this product be handled with the same precautions as used for patient samples.
- Do not use beyond the expiration date.

Preparation and Storage

- The calibrators are provided ready to use.
- Bring calibrator to room temperature (18 - 25 °C) for use.
- Always store the Calibrator Set in an upright position at 2 - 8 °C when not in use.

Stability

When stored unopened and refrigerated at 2 - 8 °C, the AIA-PACK AFP Calibrator Set is stable until the expiration date on the label. After opening, the calibrators should be used within 24 hours.

Procedure

Refer to the CALIBRATION PROCEDURE in the AIA-PACK section of this analyte application. For additional procedural instructions regarding calibration, refer to the AIA System Operator's Manual.

Calibration

1. Verify that both the calibrator lot and concentrations have been correctly entered into the analyzer software.
2. Add the appropriate amounts of each calibrator to sample cups (refer to the instrument worklist for the amount required in each sample cup).
3. Place the sample cups and the test cups on the analyzer as indicated on the worklist.
4. Tosoh recommends that all calibrators be run in triplicate.

Assignment of Values

The AIA-PACK AFP Calibrator Set contains assigned concentrations of Alpha-Fetoprotein. The assigned value is determined on a lot-to-lot basis and is designed to provide an assay calibration range of approximately 0.0 to 400.0 ng/mL of Alpha-Fetoprotein. The calibrators in this set have been standardized against WHO 1st IS 72/225 (1974).

Results

1. The mean rate for the zero calibrator should be <3.0 nM/sec.
2. Since there is a direct relationship between concentration and rate, the rates should increase as the concentration increases.
3. The replicate values should be within a 10% range.

Limitations

The AIA-PACK AFP Calibrator Set is designed solely for use with AIA-PACK assay procedures.

References

1. AIA Analyte Application Manual (AAM). Tosoh Bioscience, Inc., Grove City, OH.
2. AIA-System Operator's Manual. Tosoh Bioscience, Inc., Grove City, OH.



TOSOH CORPORATION
Bioscience Division
2-2-1 Yaesu, Chuo-ku, Tokyo 104-0028 Japan
Phone: +81 3 6636 3734
Fax: +81 3 6636 3627
Website: www.tosoh.com

Supplied by

TOSOH BIOSCIENCE, INC.
3600 Gantz Road,
Grove City, OH 43123 USA
Phone: +1 650 615 4970
Fax: +1 650 615 0415
Phone: (800) 248 6764
Fax: (800) 685 7595
Website: www.tosohbioscience.us

 European Conformity	 <i>In vitro</i> diagnostic medical device	 Consult instructions for use	 Temperature limitation
 Batch code / Lot number	 Manufacturer	 Authorized representative in the European Community	 Use by date
 Catalogue number / Part number	 Supplied by	 Sufficient for	

AIA-PACK

AFP Sample Diluting Solution

Intended Use

The AIA-PACK AFP Sample Diluting Solution is intended for IN VITRO DIAGNOSTIC USE ONLY to dilute patient samples that have concentrations of Alpha-Fetoprotein above the linear range of the assay.

Summary and Explanation

The AIA-PACK AFP Sample Diluting Solution contains human serum with no detectable concentration of Alpha-Fetoprotein. This Sample Diluting Solution is to be used only with samples that are being tested for Alpha-Fetoprotein concentrations using the ST AIA-PACK AFP assays.

Material Provided (Cat. No. 0020552)

4 x 4 mL Sample Diluting Solution

Human serum containing no detectable concentration of Alpha-Fetoprotein (0 ng AFP/mL), and 0.1% sodium azide as a preservative.

Warnings and Precautions

- The AIA-PACK AFP Sample Diluting Solution is for in vitro diagnostic use.
- These materials contain sodium azide, which may react with lead or copper plumbing to form potentially explosive metal azides. When disposing of such reagents, always flush with large volumes of water to prevent azide build-up.
- Human sera used in the preparation of this product has been tested by FDA cleared methods and found negative for the presence of HBsAg and antibody to HIV-1 and HCV. Because no test method can offer complete assurance that products derived from human blood will not transmit infectious agents, it is recommended that this product be handled with the same precautions as used for patient samples.
- Do not use beyond the expiration date.

Preparation and Storage

- The AIA-PACK AFP Sample Diluting Solution is provided ready to use.
- Always store the Sample Diluting Solution in an upright position at 2 - 8 °C when not in use.

Stability

When stored unopened and refrigerated at 2 - 8 °C, the AIA-PACK AFP Sample Diluting Solution is stable until the expiration date on the label. After opening, the Sample Diluting Solution is stable for up to 90 days when refrigerated at 2 - 8 °C.

Procedure

- Refer to the AIA System Operator's Manual for additional procedural instructions regarding sample dilution.
- If a specimen is found to contain greater than the linearity limit of 400 ng/mL, the specimen should be diluted with the Sample Diluting Solution and assayed according to the procedure in the AIA-PACK section of the analyte application.
- The AIA-900 and AIA-2000 will perform dilutions automatically if the dilution factors are entered into the software prior to assaying the diluted sample.
- The recommended dilution for serum containing greater than 400 ng/mL is 1:10 or 1:100. However, it is desirable to dilute the serum samples that contain more than 400 ng AFP/mL so that the diluted sample reads between 10 and 400 ng AFP/mL.

Results

When an auto-dilution is performed, the AIA instrument will calculate the final result.

Limitations

The AIA-PACK AFP Sample Diluting Solution is designed solely for use with AIA-PACK assay procedures.

References

1. AIA Analyte Application Manual (AAM). Tosoh Bioscience, Inc., Grove City, OH.
2. AIA-System Operator's Manual. Tosoh Bioscience, Inc., Grove City, OH.



TOSOH CORPORATION
Bioscience Division
2-2-1 Yaesu, Chuo-ku, Tokyo 104-0028 Japan
Phone: +81 3 6636 3734
Fax: +81 3 6636 3627
Website: www.tosoh.com

Supplied by

TOSOH BIOSCIENCE, INC.
3600 Gantz Road,
Grove City, OH 43123 USA
Phone: +1 650 615 4970
Fax: +1 650 615 0415
Phone: (800) 248 6764
Fax: (800) 685 7595
Website: www.tosohbioscience.us



European Conformity



In vitro diagnostic medical device



Consult instructions for use



Temperature limitation



Batch code / Lot number



Manufacturer



Authorized representative
in the European Community



Use by date



Catalogue number
/ Part number



Supplied by



Sufficient for

US: For prescription use only

AIA-PACK

AFP Calibration Verification/Linearity Test Set

Intended Use

The AIA-PACK AFP Calibration Verification/Linearity Test Set is intended for IN VITRO diagnostic use only for the verification of the calibration and linearity of the ST AIA-PACK AFP Assays.

Summary and Explanation

The AIA-PACK AFP Calibration Verification/Linearity Test Set can be used as part of a quality assurance program to assist in complying with various regulatory requirements under which an institution may operate. Specific guidelines for the acceptability of data are determined by the institution in conformance with these requirements. Good laboratory practices stipulate using different materials for calibration and quality control purposes. Calibration Verification Materials should be analyzed as unknown test samples according to the directions stated in the AIA System Operator's Manual.

Material Provided (Cat. No. 0020652)

2 x 4 mL	Sample Diluting Solution (SDS) Human serum containing no detectable concentration of AFP (0 ng AFP/mL), and 0.1% sodium azide as preservative. (Ready to Use)
2 x 2 mL	Calibration Verification Material (CVM) Human serum containing the assigned concentration of AFP (approximately 400 ng/mL, described on each vial) and 0.1% sodium azide as preservative. (Ready to Use)

Warnings and Precautions

- The AIA-PACK AFP Calibration Verification/Linearity Test Set is for in vitro diagnostic use.
- These materials contain sodium azide, which may react with lead or copper plumbing to form potentially explosive metal azides. When disposing of such materials, always flush with large volumes of water to prevent azide build-up.
- Human sera used in the preparation of this product has been tested by FDA cleared methods and found negative for the presence of HBsAg and antibody to HIV-1 and HCV. Because no test method can offer complete assurance that products derived from human blood will not transmit infectious agents, it is recommended that this product be handled with the same precautions as used for patient samples.
- Do not use beyond the expiration date.

Preparation and Storage

- The Calibration Verification Material and Sample Diluting Solution are provided ready to use.
- The materials should be at room temperature (18 - 25 °C) prior to use.
- Always store the Calibration Verification/Linearity Test Set in an upright position at 2 - 8 °C when not in use.

Stability

When stored unopened and refrigerated at 2 - 8 °C, the AIA-PACK AFP Calibration Verification/Linearity Test Set is stable until the expiration date on the label. After opening, the materials should be used within 24 hours.

Procedure

The frequency of performing calibration verification and linearity of the AIA-PACK AFP assay is established by each institution in accordance with its quality assurance program and appropriate regulatory requirements. Refer to the ASSAY PROCEDURE in the AIA-PACK section of this analyte application. For additional procedural instructions regarding auto dilution, refer to the AIA System Operator's Manual.

- Make the desired dilutions of the Calibration Verification Material with Sample Diluting Solution and mix well. A sufficient amount should be made to assay each diluted sample in triplicate.
- Program the instrument to run each of the prepared dilutions three times.
- Add the appropriate amounts of each diluted sample to sample cups (refer to the instrument worklist for the amount required in each sample cup).
- Verify the appropriate amount of test cups for the assay.
- Load sample cups on the instrument as instructed in the AIA System Operator's Manual.
- Start the assay.
- Record the values obtained on the appropriate forms and calculate the desired statistics.

Assignment of Values

The Calibration Verification Material contains an assigned concentration of Alpha-Fetoprotein. The assigned value is determined on a lot-to-lot basis and is designed to approximate the upper linear range of the assay. The Calibration Verification Material has been standardized against WHO 1st IS 72/225 (1974).

Results

- Make a determination of the acceptability of the data according to the specific guidelines established by your institution which satisfy the regulatory requirements under which your institution operates. If results do not meet your specifications, please initiate corrective action, as appropriate, or contact Tosoh Bioscience for assistance.
- Retain test records in the laboratory in the designated location for future reference.

Limitations

The AIA-PACK AFP Calibration Verification/Linearity Test Set is designed solely for use with AIA-PACK assay procedures. The AIA-PACK AFP Calibration Verification/Linearity Test Set will only verify linearity between the established sensitivity of the assay and the value of analyte level listed on the Calibration Verification Material label.

References

1. AIA Analyte Application Manual (AAM). Tosoh Bioscience, Inc., Grove City, OH.
2. AIA-System Operator's Manual. Tosoh Bioscience, Inc., Grove City, OH.



TOSOH CORPORATION
Bioscience Division
2-2-1 Yaesu, Chuo-ku, Tokyo 104-0028 Japan
Phone: +81 3 6636 3734
Fax: +81 3 6636 3627
Website: www.tosoh.com

Supplied by

TOSOH BIOSCIENCE, INC.
3600 Gantz Road,
Grove City, OH 43123 USA
Phone: +1 650 615 4970
Fax: +1 650 615 0415
Phone: (800) 248 6764
Fax: (800) 685 7595
Website: www.tosohbioscience.us

 European Conformity	 In vitro diagnostic medical device	 Consult instructions for use	 Temperature limitation
 Batch code / Lot number	 Manufacturer	 Authorized representative in the European Community	 Use by date
 Catalogue number / Part number	 Supplied by	 Sufficient for	