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
Figure 2. Impurities formed as a result of HC mispairing (light blue dotted box), LC mispairing (light red dotted box), fragmentation (light green dotted box) and aggregation (light purple dotted box) are illustrated here for a representative desired asymmetric bsAb, with the proposed strategies to remove them shown in the corresponding dark colored non-dotted boxes. Not shown here: (1) for a Fab × scFv bsAb, differential KappaSelect affinity and CH1-based chromatographic methods have proven to be useful [59, 60]; (2) diabody-IgG mispaired products can be separated from scFv-IgG bsAb targets through the use of an alkaline pH near the pI of the bsAb using cation exchange chromatography [71]; (3) hydrophobic and mixed-mode resins have been proposed to provide good separation for certain homodimer mispaired products from their bsAb targets [77, 78]. Examples listed here are not exhaustive.

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
文献 1 : S. W. Chen et al., Current trends and challenges in the downstream purification of bispecific antibodies, Antibody Therapeutics, 2021, 4, 73-88, CC BY,

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1. S. W. Chen et al., Current trends and challenges in the downstream purification of bispecific antibodies, *Antibody Therapeutics*, 2021, 4, 73-88, <https://doi.org/10.1093/abt/tbab007>
2. X. Chen et al., Protein L chromatography, A useful tool for monitoring/separating homodimers during purification of IgG-like asymmetric bispecific antibodies, *Protein Expression and Purification*, 175 (2020) 105711, <https://doi.org/10.1016/j.pep.2020.105711>
3. C. Chen et al., Controlled conductivity at low pH in Protein L chromatography enables separation of bispecific and other antibody formats by their binding valency, *MABS*, 2019, 11, 632-638, <https://doi.org/10.1080/19420862.2019.1583996>
4. U. Brinkmann et al., The making of bispecific antibodies, *MABS* 2017, 9, 182-212, BY NC ND, <https://doi.org/10.1080/19420862.2016.1268307>
5. S. M. Chamow et al., Capture of CH1-containing bispecific antibodies, *BioProcess International*, 18(5), May 2020, <https://bioprocessintl.com/downstream-processing/chromatography/capture-of-ch1-containing-bispecific-antibodies-evaluating-an-alternative-to-protein-a/>
6. C. Andrade et al., An integrated approach to aggregate control for therapeutic bispecific antibodies using an improved three column Mab platform-like purification process, *Biotechnol. Progress*, 2019, 35, 1-13, BY NC ND, <https://doi.org/10.1002/btpr.2720>
7. Process analytics and intermediate purification of bispecific antibodies with a non-affinity platform, *Tosh Application Note*, a19p09a, <https://www.separations.eu.tosohbioscience.com/OpenPDF.aspx?path=~/File%20Library/TBG/Products%20Download/Application%20Note/a19p09a.pdf>
8. L. K. Kimerer et al., Chromatographic behavior of bivalent bispecific antibodies on hydrophobic interaction chromatography columns, *J. Chromatogr., A*, 1617 (2020) 460836, <https://doi.org/10.1016/j.chroma.2019.460836>
9. Y. Li et al., A roadmap for IgG-like bispecific antibody purification, *Approaches to the Purification, Analysis and Characterization of Antibody-Based Therapeutics*, 2021, 167-179, <https://doi.org/10.1016/B978-0-08-103019-6.00008-4>
10. Y. Wan et al., Application of pH-salt dual gradient elution in purifying a WuXiBody-based bispecific antibody by MMC ImPres mixed-mode chromatography, *Protein Expression and Purification*, 181 (2021) 105822, <https://doi.org/10.1016/j.pep.2021.105822>
11. S. W. Chen et al., Excellent removal of knob-into-hole bispecific antibody byproducts and impurities in a single-capture chromatography, *Bioresources and Bioprocessing*, 2022, 9:72, CC BY, <https://doi.org/10.1186/s40643-022-00562-y>
12. J. A. Waller et al., Ceramic hydroxyapatite chromatography plays a critical role in bispecific antibody purification process for impurity removal, *Antibody Therapeutics*, 2023, 6, 30-37, <https://doi.org/10.1093/abt/tbac030>
13. 東ソー、テクニカルノート TOYO No. 2 「抗体精製用充填剤は、これで決まり！！(2)」

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