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

S. Kurth et al., Separation of full and empty adeno-associated virus capsids by anion-exchange chromatography using choline-type salts, *Analytical Biochemistry* 686 (2024) 115421, CC BY-NC-ND, <https://doi.org/10.1016/j.ab.2023.115421>

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カラム ; TSKgel Q-STAT (7 μm, 4.6 mm I.D. x 10 cm)
 溶離液 ; Buffer A: 20 mmol/L Tris-HCl (pH 9.0)
 ; Buffer B: Buffer A containing 1 mol/L choline chloride
 ; 20 min linear gradient from Buffer A/B (90/10) to (65/35)
 流速 ; 1.0 mL/min
 検出 ; Fluorescent, excitation at 280 nm, emission at 350 nm
 温度 ; 30 °C
 試料 ; purified AAV5 (5 x 10¹² vc/mL) of empty/full (50/50) was repeatedly diluted to 1/4 for sample preparation

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
1. 東ソー、テクニカルノート TSKgel No. 15、ウイルスベクターの分離、分析はこれで決まり！
2. 東ソー、テクニカルノート TSKgel No. 20、アデノ随伴ウイルス (AAV) の分離、分析はこれで決まり！
3. R. Rieser et al., Comparison of different liquid chromatography-based purification strategies for adeno-associated virus vectors, *Pharmaceutics*, 13 (2021) 748, <https://doi.org/10.3390/pharmaceutics13050748>
4. K. Richter et al., Purity and DNA content of AAV capsids assessed by analytical ultracentrifugation and orthogonal biophysical techniques, *Eur. J. Pharm. Biopharm.*, 189 (2023) 68-83, <https://doi.org/10.1016/j.ejpb.2023.05.011>
5. Z. Wu et al., Development of a two-dimensional liquid chromatography-mass spectrometry platform for simultaneous multi-attribute characterization of adeno-associated viruses, *Anal. Chem.* 24 (2022) 3219-3226, <https://doi.org/10.1021/acs.analchem.1c04873>
6. M. K. Aebischer et al., Anion-exchange chromatography at the service of gene therapy: baseline separation of full/empty adeno-associated virus capsids by screening of conditions and gradient elution mode, *J. Mol. Sci.*, 23 (2022) 12332, <https://doi.org/10.3390/jms232012332>
7. P. L. R. H. Joshi et al., Development of a scalable and robust AEX method for enriched rAAV preparations in genome-containing VCs of serotypes 5, 6, 8, and 9, *Molecular Therapy: Methods & Clinical Development*, 21 (2021) 341, <https://doi.org/10.1016/j.omtm.2021.03.016>
8. S. L. Khatwani et al., Anion-exchange HPLC assay for separation and quantification of empty and full capsids in multiple adeno-associated virus serotypes, *Molecular Therapy: Methods & Clinical Development* 21 (2021) 548-558, <https://doi.org/10.1016/j.omtm.2021.04.003>
9. C. Wang et al., Developing an Anion Exchange Chromatography Assay for Determining Empty and Full Capsid Contents in AAV6.2, *Molecular Therapy: Methods & Clinical Development*, 15 (2019) 257-263, <https://doi.org/10.1016/j.omtm.2019.09.006>
10. M. Mietzsch et al., Production and characterization of an AAV1-VP3-only capsid: An analytical benchmark standard, *Molecular Therapy: Methods & Clinical Development*, 29 (2023) 460-472, <https://doi.org/10.1016/j.omtm.2023.05.002>
11. J. Kervinen, Anion-exchange chromatography separation of full and empty capsids using non-toxic choline-type salts, BPN Conference, Australia, 2023

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