




TOSOH

**Quality Immunoassays  
for Quality of Life**

**Clinical Testing Solutions with  
Consistency, Stability and Harmony**



### **Unit Dose Test Cup Reagent System**

- Dry reagent format
- 90-day calibration stability for most assays
- Single, unitized cups - No pre-mixing or measuring
- Interchangeable reagents between Tosoh AIA® Systems
- Barcoded for identification/inventory management
- 10 minute incubation time for most assays
- Free of biotinylated antibodies by design
- Extensive test menu includes assays for cardiac, thyroid, reproductive, tumor, metabolic, diabetes, kidney, and anemia markers, Vitamin D and more



**TOSOH BIOSCIENCE**

[www.tosohbioscience.us](http://www.tosohbioscience.us)

# Tosoh AIA Test Menu

10-minute incubation unless otherwise noted.

## Tumor Markers

AFP  
CEA  
PA  
27.29  
CA125  
CA19-9

## Diabetes

C-Peptide II  
IRI

## Thyroid

TSH  
TT3  
T4  
FT3  
FT4  
T-U  
TSH 3rd-Gen ♦  
TPOAb ♦  
TgAb ♦

## Cardiac Markers

BNP\*  
CK-MB  
cTnI 2nd-Gen  
Myoglobin

## Anemia

FER  
B12 ♦  
FOLATE ♦  
RBC FOLATE ♦

## Reproductive

βHCG  
HCG  
DHEA-S  
E2  
FSH  
LH II  
PRL  
PROG II  
PROG III  
SHBG  
Testosterone

## Kidney Markers

BMG  
Cystatin C  
intact PTH

## Metabolic

CORT  
HGH

## Additional

ACTH  
Homocysteine  
IgE II  
PAP  
25-OH Vitamin D ♦

♦ > 10 minute time for pretreatment and/or incubation. Not available on the AIA-360.

\* For Emergency Department use

## One reagent for all Tosoh AIA Systems

The AIA-PACK® Test Cup format works with every Tosoh AIA Automated Immunoassay System. Laboratories with more than one Tosoh system can use the same reagents interchangeably between instruments. Transition between systems is seamless, ensuring consistent results and efficient, economical operation for the laboratory.



**AIA-360**



**AIA-900**  
Benchtop



**AIA-900**  
Loader Model



**AIA-900**  
with 9 Tray Sorter  
(Also available in 19 tray  
sorter option)



**AIA-2000**

Not all immunoassays are the same. Using the right immunoassay gives you confidence in the results you deliver.

## Choose wisely. Choose Tosoh.



# TOSOH BIOSCIENCE

TOSOH

Tosoh Bioscience, Inc.  
3600 Gantz Road  
Grove City, OH 43123  
Tel: (800) 248-6764  
Fax: (650) 615-0415

[www.tosohbioscience.us](http://www.tosohbioscience.us)

*"Tosoh Bioscience", "AIA-PACK" and "AIA" are the registered trademarks of Tosoh Corporation*