




Title: Assigned Values for Hemoglobin A1c Control Set Part #0021974 IFU, U.S. Only

Table Revised: April, 2026

Hemoglobin A1c Control Set Part #0021974*					
LOT	Level 1		Level 2		 Use-by date
	IFCC	NGSP	IFCC	NGSP	
AB5060	34 ± 3 mmol/mol	5.3 ± 0.3%	87 ± 5 mmol/mol	10.1 ± 0.5%	2027-11-30 (YYYY-MM-DD)

*For use in the U.S. with Tosoh Automated Glycohemoglobin Analyzer HLC-723GR01