

# AIA-PACK

## Sample Treatment Cup

---

### Intended Use

The AIA-PACK Sample Treatment Cup is intended for IN VITRO DIAGNOSTIC USE ONLY. The AIA-PACK Sample Treatment Cup can be used for automatic sample dilutions and automatic sample pretreatment with any AIA-PACK Assay on the Tosoh AIA-900 and AIA-2000 Immunoassay Analyzers.

### Material Provided (Cat. No. 0020971)

10 trays x 20 cups - (Each tray contains 20 sample treatment cups)

### Warnings and Precautions

The AIA-PACK Sample Treatment Cup is intended for in vitro diagnostic use only.

### Preparation and Storage

- The AIA-PACK Sample Treatment Cup is provided ready to use.
- Store at room temperature (18 ° - 25 ° C).

### Stability

- The AIA-PACK Sample Treatment Cup is stable until the expiration date on the label.

### References

1. AIA Analyte Application Manual (AAM). Tosoh Bioscience, Inc., Grove City, OH.
2. AIA System Operator's Manual. Tosoh Bioscience, Inc., Grove City, OH.

# AIA-PACK

## Sample Treatment Cup

---



**TOSOH CORPORATION**  
**Bioscience Division**  
2-2-1 Yaesu, Chuo-ku, Tokyo 104-0028 Japan  
Phone: +81 3 6636 3734  
Fax: +81 3 6636 3627  
Website: [www.tosoh.com](http://www.tosoh.com)

Supplied by

**TOSOH BIOSCIENCE, INC.**  
3600 Gantz Road,  
Grove City, OH 43123 USA  
Phone: +1 650 615 4970  
Fax: +1 650 615 0415  
Phone: (800) 248 6764  
Fax: (800) 685 7595  
Website: [www.tosohbioscience.us](http://www.tosohbioscience.us)

|                                                                                                                        |                                                                                                                                 |                                                                                                                                               |                                                                                                               |
|------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------|
| <br>European Conformity               | <br><i>In vitro</i> diagnostic medical device | <br>Consult instructions for use                           | <br>Temperature limitation |
| <br>Batch code / Lot number           | <br>Manufacturer                             | <br>Authorized representative<br>in the European Community | <br>Use by date            |
| <br>Catalogue number<br>/ Part number | <div>Supplied by</div>                                                                                                          | <br>Sufficient for                                         | <div>US: For prescription use only</div>                                                                      |