



AIA-360

Best Tests for Pets



AIA-360

Accurate Immunoassay Results with User-friendly Operations

- Fast turnaround ~ 20 minutes; throughput of 36 tests/hour
- Carousel holds 25 samples and 25 test cups at one time
- Compact - 16"w x 16"d x 21"h / 61 lbs
- Random access
- Once samples and cups are loaded, just press "Start" and walk away
- Extensive test menu includes Progesterone, T4, Insulin

TOSOH BIOSCIENCE

www.tosohbioscience.us

Accurate Results with User-Friendly Operations

AIA-360 Provides Simplified and Efficient Immunoassay Testing




Place Test Cups



Place Samples



Press Start



Unit Dose Test Cup Reagents

- Single, unitized test cups require no pre-mixing, no pre-measuring to minimize preparation time and waste
- Dry reagent format ensures 90 day calibration stability for most assays
- Bar-coded for easy test identification and inventory management
- All Tosoh AIA analyzers utilize the same reagents



AIA-360
36 Tests Per Hour



AIA-900 with 9 Tray Sorter
90 Tests Per Hour
*Also available with 19 Tray Sorter
and as Loader Model*



AIA-2000
200 Tests Per Hour



TOSOH BIOSCIENCE

"Tosoh Bioscience" and "AIA" are the registered trademarks of Tosoh Corporation

Tosoh Bioscience, Inc.
3600 Gantz road
Grove City, OH 43123
Tel: (800) 248 6764
Fax: (650) 615 0415

www.tosohbioscience.us