

AIA-PACK

Detector Standardization Test Cup

Intended Use

The AIA-PACK Detector Standardization Test Cup is intended for IN VITRO DIAGNOSTIC USE ONLY. The AIA-PACK Detector Standardization Test Cup is used to measure the substrate background and calibrate the detectors on the Tosoh AIA System Analyzers

Material Provided (Cat. No. 0020970)

10 trays x 20 cups (Each tray contains 20 detector standardization cups)

Warnings and Precautions

- The AIA-PACK Detector Standardization Test Cup is for in vitro diagnostic use only.

Preparation and Storage

- The AIA-PACK Detector Standardization Test Cup is provided ready to use.
- Store at room temperature (18 ° - 25 ° C).

Stability

- The AIA-PACK Detector Standardization Test Cup is stable until the expiration date on the label.

AIA-PACK

Detector Standardization Test Cup

Procedure

For additional procedural instructions regarding detector standardization and substrate background measurement, refer to the AIA System Operator's Manual.

References

1. AIA Analyte Application Manual (AAM). Tosoh Bioscience, Inc., Grove City, OH.
2. AIA System Operator's Manual. Tosoh Bioscience, Inc., Grove City, OH.

	TOSOH CORPORATION Bioscience Division 2-2-1 Yaesu, Chuo-ku, Tokyo 104-0028 Japan Phone: +81 3 6636 3734 Fax: +81 3 6636 3627 Website: www.tosoh.com
Supplied by	TOSOH BIOSCIENCE, INC. 3600 Gantz Road, Grove City, OH 43123 USA Phone: +1 650 615 4970 Fax: +1 650 615 0415 Phone: (800) 248 6764 Fax: (800) 685 7595 Website: www.tosohbioscience.us

 European Conformity	 In vitro diagnostic medical device	 Consult instructions for use	 Temperature limitation
 Batch code / Lot number	 Manufacturer	 Authorized representative in the European Community	 Use by date
 Catalogue number / Part number	Supplied by	 Sufficient for	US: For prescription use only