




Title: Assigned Values for HbA1c Calibrator Set (S) Part #0023502 IFU, U.S. Only

Table Revised: April, 2026

HbA1c Calibrator Set (S) Part # 0023502*					
LOT	HbA1c Calibrator (1)		HbA1c Calibrator (2)		 Use-by date
	IFCC	NGSP	IFCC	NGSP	
ZB5030	41.1 mmol/mol	5.91 %	97.6 mmol/mol	11.08 %	2027-11-30 (YYYY-MM-DD)

\*For use in the U.S. with Tosoh Automated Glycohemoglobin Analyzer HLC-723GR01