



IVF Solutions

E2 Results in 18 to 20 Minutes



Why Wait?

Tosoh Immunoassay Analyzers Deliver Accurate Results Quickly

- Fast** First Result in 18 to 20 Minutes
- Efficient** Simple Automated Testing
- Accurate** Reliable Results
- Consistent** Common Unit Dose Test Cup Reagent
- Cost Saving** 90 Day Calibration Stability for all Fertility Assays
- Comprehensive** Full Immunoassay Menu Available, Including a Full Menu of Reproductive Hormone Assays
- Biotin Free Immunoassay Design** Tosoh's Immunoassays are free of biotinylated antibodies



Tosoh's fast results allow review of ultra-sound and E2 on the same visit.

TOSOH BIOSCIENCE

Unit Dose Test Cup Reagent Technology

- Single, unitized cups use a dry reagent format that requires no pre-mixing, no pre-measuring, no on-board refrigeration and no reagent waste
- Dry reagent format provides calibration stability of up to 90 days for most assays
- Reagents have a one year shelf life from date of manufacture
- The AIA® test cups and trays are labeled with code and lot number for automated scheduling and inventory management on AIA-2000 and AIA-900 sorter models
- ST AIA-PACK® reagents provide a fast turnaround time due to a 10 minute incubation protocol
- Extensive test menu provides a wide range of immunoassay tests including thyroid, cardiac markers, tumor markers, reproductive hormones, anemia markers, metabolic markers, diabetes, kidney markers and more



Available for Labs and Workloads of All Sizes



AIA-360
36 Tests Per Hour



AIA-900 Benchtop
90 Tests Per Hour



AIA-900 Loader Model
90 Tests Per Hour



AIA-900 with 9 Tray Sorter
90 Tests Per Hour
(Also available in 19 tray sorter option)



AIA-2000
200 Tests Per Hour

Choose Wisely. Choose Tosoh.



TOSOH BIOSCIENCE

TOSOH

Tosoh Bioscience, Inc.
3600 Gantz Rd
Grove City, OH 43123
Tel: (800) 248-6764
Fax: (650) 615-0415
www.tosohbioscience.us

"AIA-PACK" and "AIA" are registered trademarks of Tosoh Corporation