



TOSOH

Vitamin D

AIA-PACK® Series

A New Day
for Vitamin D



Fortify Your Lab with Vitamin D Testing

Vitamin D measurement is intended as an aid in the determination of Vitamin D sufficiency.

The antibodies have equimolar affinity for Vitamin D2 and D3 and the assay measures Total Vitamin D. Dry reagent format ensures 90-day calibration stability for minimal waste and cost effective testing.

TOSOH BIOSCIENCE

www.tosohbioscience.us

Universal Reagent

Tosoh's AIA-PACK test cup format works with every Tosoh automated immunoassay system allowing for a seamless transition from one system to the other, ensuring consistent results in an efficient and economical process.

Save Time and Money

Tosoh's AIA-PACK test cups are single, unitized cups that use a dry reagent format that ensures calibration stability of up to 90 days.

Minimal Contamination Risk, More Traceability

Because there is no transfer of reagents the risk of contamination is eliminated. The unique bichromatic fluorescence kinetic measurement ensures a high analytical and functional sensitivity for all assays. AIA-PACK test cups and trays are labelled with the assay code and lot number for automated scheduling and inventory.

Fortify Your Lab with Vitamin D Testing

- The Tosoh ST AIA-PACK 25-OH Vitamin D correlates well with other immunoassay methods and measures 25-OH Vitamin D2 and D3 in equimolar proportions
- The ST AIA-PACK 25-OH Vitamin D assay is aligned to the reference measurement procedure used in the Vitamin D standardization program VDSP for Vitamin D measurement
- 90 day calibration stability makes it feasible to run this test even with low volumes
- Available on a broad range of analyzers – AIA-900 and AIA-2000
- Utilizes Tosoh's Unit Dose Test Cup reagent technology
- Assay range of 4 to 120 ng/mL
- Sample type - Serum, Na heparinized plasma or EDTA plasma

Available for Labs and Workloads of All Sizes



Part #	Product Name	Description
0025234	ST AIA-PACK 25-OH Vitamin D	20 cups x 5 trays (100 tests)
0025334	ST AIA-PACK 25-OH Vitamin D CALIBRATOR SET	Six levels, 2 x 6 x 1 mL
0025734	ST AIA-PACK 25-OH Vitamin D PRETREATMENT SET	PR-1-1 x 32 mL, 6 Reagent vials (empty) PR-2-6 x 5 mL
0025434	AIA-PACK 25-OH Vitamin D CONTROL SET	Two Levels, 2 x 2 x 1 mL

Tosoh products are for Prescription Use Only as In-Vitro Diagnostics



Tosoh Bioscience, Inc.
 3600 Gantz Road
 Grove City, OH 43123
 Tel: (800) 248-6764
 Fax: (650) 615-0415
www.tosohbioscience.us

"Tosoh Bioscience", "AIA-PACK" and "AIA" are the registered trademarks of Tosoh Corporation